



**Soleoline**

# SOLEO LINE

## SOLEO LINE

## Ultrasound, electro- and combination therapies

**Soleo Sono, Soleo Galva und Soleo SonoStim**  
These three therapeutic devices from the Soleo-line series are innovative, powerful and very easy to use thanks to their small size. Ideal for mobile use. Whether in your practice, during house calls, or at sporting events, you can use them wherever your patients are.

### **Extremely easy to use**

- A large, colour display with one-touch operation
- Easy programme selection based on treatment goal
- Easy to use thanks to intuitive menu navigation



**Expertise for any medical problem**

Clearly presented, recommended treatments for various therapies available at all times.

**Tailored to your practice**

- The programme's customised start settings save you precious treatment time.
- Change pre-set programmes and save them to fit your patients.
- Your custom settings are saved to an SD card.

The comprehensive design of the new *Soleoline* generation combines aesthetics with functionality, ergonomics and durability.

**Simply beautiful.**

A sleek case and an attractive, capacitive glass touchscreen ensure that *Soleo Sono*, *Soleo Galva* und *Soleo SonoStim* really stand out in your practice.

**Beautifully simple.**

The user interface of the new *Soleoline* generation has been redesigned to be clear and self-explanatory. It's even fun to use with every application.

Intelligent design –  
sets new standards



# SOLEO SONO

Soleo Sono combines the classic possibilities of modern ultrasound therapy with the new SonoSwing process.

#### **Helps wherever the pain is**

The ergonomic, multiple frequency transducer can give off two frequencies: 0.8 MHz (for penetration) and 2.4 MHz (to treat the surface).

#### **Unique**

The patented SonoSwing process allows you to adjust the percentage of both frequencies. Due to the adjustable mix of both frequencies you can choose the depth of the optimal heat development.

#### **Time-saving and safe**

- Select treatment based on symptoms
- Recommended treatments for 50+ symptoms
- Up to 120 favourite and customised programmes can be stored

## The new dimension of ultrasound therapy

# For the full spectrum of modern electrotherapy

With a wide variety of current types, special sports and muscle rehabilitation programmes Soleo Galva offers a comprehensive range of treatments.

## Two channels

The number of different varieties in the 2-channel mode increases your range of therapies.

Twinstim: 1 programme / 2 channels  
Duostim: 2 programmes / 2 channels

## Sports therapy

Special additional current types and programmes allow you to expand the different types of sports therapy you offer.

## VacoS: Electrotherapy plus suction massage

VacoS allows you to apply electrodes with ease. Your patients are sure to enjoy the therapeutic effects of a suction massage during electrotherapy.



# SOLEO SONOSTIM

## SOLEO SONOSTIM

## Ultrasound and electrotherapy all in one

Soleo *SonoStim* does everything that Soleo *Sono* and Soleo *Galva* do, and it also allows you to combine ultrasound therapy with electrotherapy. With its innovative technology, more than 300 indications, this mobile and high-tech device has a huge range of applications.

Zimmer  
MedizinSystem



### **Simultaneous therapy**

The simultaneous delivery of ultrasound and electrotherapy intensifies the therapeutic effects of each.

Ultrasound and electrotherapy combined  
- Really easy to use

Due to an adjustable mix of both frequencies you determine the optimal depth of heat Generation.

# High-Performance Therapy

## Technical Data



### Ultrasound therapy version

Ultrasound head

Frequencies

Maximum power emission

User guide

### *Sono*

1 socket for one 5 cm<sup>2</sup> multi-frequency watertight head, (1 cm<sup>2</sup> as an optional extra)

0.8 MHz: deep effect / 2.4 MHz: superficial effect

3 W/cm<sup>2</sup> for the 5 cm<sup>2</sup> head

More than 50 therapy recommendation

### Electrostimulation version

Channels

User guide

### *Galva*

2 independent channels (programmes, intensities, parameters)

More than 120 therapy recommendation

### Combined version

General features

User guide

### *SonoStim = Sono + Galva*

Same features as *Sono* and *Galva*

More than 300 therapy recommendation

### *Soleoline*

SD memory card

Dimensions

### **For all versions**

120 favourites entries

326 mm x 134 mm x 243 mm / 2,1 kg

### Vaco equipment

Operating modes

Control

Dimensions

### **VacoS for *Galva* or *SonoStim* version**

Continuous and pulsed

via the touch screen

322 mm / 234 mm / 130 cm / Weight 3 kg

Zimmer MedizinSysteme GmbH  
 Junkersstraße 9  
 89231 Neu-Ulm, Germany  
 Tel. +49 7 31. 97 61-291  
 Fax +49 7 31. 97 61-299  
 export@zimmer.de  
 www.zimmer.de

