

Lokomat®

ORTHOPEDIC MATERIAL STORAGE

The Lokomat portfolio provides two orthopedic storage solutions, to store your orthopedic material, which allow ergonomic handling and faster training preparation. The orthopedic material storage solutions are optional accessories to the LokomatNanos and the LokomatPro.



All content is subject to change without notice. © Hocoma AG, CH-Lokomat-OrthopedicStorage 171107 EN

Trolley for the safe and easy storage of orthopedic material as well as up to four orthoses (two pairs)

Measures:

88 cm (L) x 78 cm (W) x 123.4 cm (H) 34.65 in (L) x 30.71 in (W) x 48.58 in (H)





Wall mounting system for the safe and easy storage of the most frequently used orthopedic material (cuffs, harnesses) as well as up to two additional orthoses (one pair)

Measures:

70 cm (W) x 110.5 cm (H) x 31.9 cm (D) 27.56 in (W) x 43.5 in (H) x 12.56 in (D)

Key features

- 1. Safe and easy storage
- 2. Easier training preparation
- 3. Faster patient set up
- 4. More ergonomical handling of orthopedic material and orthoses

For more information on prices and financing options,

please do not hesitate to contact us.

CONTACT US

Hocoma AG, Industriestrasse 4 8604 Volketswil, Switzerland Phone: +41 43 444 22 00 E-mail: info@hocoma.com www.hocoma.com



PRODUCT DISCLAIMER

All Hocoma products are medical devices and must be used in strict adherence to the User Manual; failure to do so may result in serious personal injury. It is strongly recommended that you regularly consult Hocoma's website (www.hocoma.com/legalnotes) for the latest available information. Please contact Hocoma should you have any questions. Use only under the supervision of qualified medical personnel. However, certain Hocoma products are marketed for home use and must be strictly used according to the recommendations of your medical care provider who is knowledgeable about your specific needs. Consult the User Manual and Hocoma's website (www.hocoma.com/legalnotes) for appropriate product designation. Failure to obtain and follow the recommendations of your medical care provider may result in serious personal injury. This information provides details about medical products which may not be available in all countries and may not have received approval or market clearance by all governmental regulatory bodies throughout the world. Nothing herein should be construed as a solicitation or promotion of any product or of an indication of any specific use for any product which is ot authorised by the laws and regulations of the country where the reader of this information resides.