

Six ways Hi-Pro 2 boosts fitting efficiency

- Universal programming interface for all hearing instrument brands
- USB connectivity always ready
- Faster than Hi-Pro predecessors
- No training required
- Sleek, space-saving design
- Backward compatible works with any hearing instrument that has been programmed by a Hi-Pro device, regardless of version



natus

Your connection to greater fitting efficiency



Streamline your fitting processes

Hi-Pro 2 gives you an easy way to streamline your fitting processes. Developed in collaboration with hearing aid manufacturers, Hi-Pro 2 is the latest in the Hi-Pro series — long recognized as the industry standard for hearing instrument programming. Hi-Pro 2 has a universal hardware interface to program all brands and types of hearing aids — making hearing instrument programming fast, easy and reliable.

Save time

Hi-Pro 2 saves you time in the fitting process. Programming is faster than its predecessor. The high-speed mode enables fitting modules to communicate with the hearing instrument up to five times faster than the previous Hi-Pro USB version.

And unlike battery-driven programmers, the Hi-Pro 2 has plugand-play USB connectivity, so it's ready — and working — when you need it. No training required.

Boost efficiency

One programming device covers your programming needs. Hi-Pro 2 programs all major hearing instruments. Sleek and discrete, Hi-Pro 2 comes as a space-saving stand-alone device or as an integrated part of the audiometer, Aurical® Aud. Hi-Pro 2 is also "backward" compatible, so it works with any hearing instrument that has been previously programmed by a Hi-Pro device, regardless of version.

TECHNICAL SPECIFICATIONS

PC Interface:

- Communication: USB 2.0 full speed (USB 1.1 compatible)
- USB Connector Type: "Type B" connector (on the Hi-Pro 2 unit)

Power Supply:

• The Hi-Pro 2 unit is powered from the PC USB port

Operating systems:

- Windows XP Pro SP3, Windows 7 (32 and 64 bit),
- Windows 8 (32 and 64 bit), Windows 10



Choose Hi-Pro 2 for universal programming and efficient fitting. Contact your local Natus representative today or visit hearing-balance.natus.com/hipro

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. P/N 7-26-99000-EN/05.

