

Therapy tables

1. Configuration Guide	page 02
2. Accessories	page 10
2.1 Cushions for posture	page 10
2.2 Electrical accessories	page 12
2.3 For protection and hygiene	page 13
2.4 Padded supports and elements	page 13
2.5 Specific accessories for Test Series tables	page 14
3. For examination, physiotherapy and manipulations	page 15
3.1 VISIT tables range	page 15
3.2 THER tables range	page 23
3.3 UNIX tables range	page 33
3.4 SINTHESI tables range	page 37
3.5 BOBATH tables range	page 45
4. Multi-purpose tables and for specific treatments	page 53
4.1 TEST tables range	page 53
5. Massage tables and for aesthetics	page 61
5.1 ALUX tables range	page 61
6. Basic tables	page 67
7. Other furniture	page 77


Therapy tables Configuration Guide

PRODUCT CODES EXPLAINED

The product code for Chinesport height-adjustable treatment tables is an “open code” that varies depending on the choice of various options available for each model. Each code position has a particular meaning, which corresponds with one or more features of the table.

It is also worth remembering that our Research & Development department is continually developing new options and accessories to achieve table configurations that increasingly meet the professionals’ individual requirements.

POSITIONS 1-2



Details on our height-adjustable therapy tables

With an experience spanning over several decades, Chinesport offers a range of height-adjustable treatment tables for various uses and with different structural characteristics. Our offering covers various options including tables for examination, transfer, physical therapy, manipulation, specific treatment and massage. The frame is always sturdy, stable, with all rotary joints such as levers and hinges running on self-lubricating bushings to preserve them from wear and stop them creaking during movement. Finally, the base is open at the sides to allow for fitting hoisting equipment or medical equipment trolleys with overall dimensions up to 18 cm.



L V

VISIT tables range

This is the classic range of examination and treatment tables consisting mainly of two sections, with a 68 cm-long head section that can be used as a backrest. The backrest can have positive and negative tilt adjustment: up to $+70^{\circ}/-20^{\circ}$. Height is adjusted via one of the most reliable motors available on the market. A hydraulic version is also available as an alternative. The frame is unique for its simplicity and harmonious design, which is protected by registered Community design. Some of our tables are also available with larger wheels so they can be used to transport patients. This configuration can be completed by fitting a set of safety siderails.



L H

THER tables range

A major difference from the VISIT range is given by the head section length which is only 52 cm. For this reason, it cannot be used as a backrest, but it can be fitted with armrests to provide greater comfort for the patient when in prone position. These tables were developed for physiotherapy, osteopathy or for the so-called “Global Postural Re-education” method, and they can be configured with up to six sections to achieve different postural positions for the patient. There is also a model that can be arranged in a seated position for particular use in gynecology. Another difference from the VISIT range is that the head section has a greater tilt adjustment, i.e. positive tilt up to $+85^{\circ}$ and negative tilt up to -35° , plus the nose/mouth breathing hole is standard.



L C

UNIX tables range

This range is an alternative to the classic ranges called VISIT and THER used for examinations, treatments and physiotherapy. In particular, the distinguishing feature in these tables is that during height adjustment the individual sections remain completely motionless. This makes them particularly suitable for smaller environments or to make the best use of available space, as well as to maintain the position chosen for exercises once the pulley therapy system has been set up. You can opt for a simplified base frame and configure the head section with a backrest that is 68 cm long, instead of a smaller section, 52 cm long, that can come with armrests. The second type of head section can have positive and negative tilt adjustment: up to $+85^{\circ}/-35^{\circ}$.



L S

SINTHESI tables range

Our Sinthesi couches allow therapists to effortlessly place patients in countless myofascial stretching postures that can be easily maintained for longer periods of time and without strain. In this relaxed position, various techniques can be applied to facilitate recovery in myofascial elongation, increase joint movement range and achieve an antalgic posture as a relief from acute pain. Also, our Sinthesi tables provide the ideal starting position for exercises that stimulate central activation and for isometric muscular work.



L B

BOBATH tables range

Large physical therapy tables that are up to two meters long and have a surface area measuring four square meters. The size and stability characteristics mean they are suitable for neurological treatments, as patients can be guided through a wide range of exercises and sensory experiences, as well as postural adjustments. This range also includes the models used in the unique C.I.R. - Continued Intensive Rehabilitation method. Finally, a special version also allows for use with static supine exercises without having to move a patient to a rehabilitation gym.



L E

TEST tables range



This is a range of multi-functional tables for examination and for specific treatments. In particular, they can be used in podiatry, gynecology, urology, cardiology, geriatrics, beauty, phlebotomy or even to transport patients within care facilities. The table height can be adjusted electrically or via hydraulic piston, while the moving sections are servo-assisted by gas spring. In some versions, the seating section is tilt-adjustable via another motor. As a result, these tables can easily change into a seated, relaxation, Trendelenburg and anti-Trendelenburg position even with the patient lying on the table. Finally, specialized accessories are available to suit individual uses.



L A

ALUX tables range

With this range, Chinesport is branching out into beauty, spa and massage tables for the first time, and introduces some additional optional features such as the sections heated at three different temperatures and the synthetic leather upholstery that covers the sides to produce a softer padding effect on the contact surface. Among other noteworthy options, there is a possibility to choose a led lights configuration for combining chromotherapy effects, and even when there is only one option of section width or padding density/thickness clients can request personalized features and options. These tables have two motors, so they can achieve the Trendelenburg and anti-Trendelenburg position via a dual telescopic column system.



POSITIONS 3-4-5

Number and type of sections

This part of the code identifies the type of section. The main difference between models has to do with the number of sections and their dimensions, as well as with the posture that patients can put themselves in when using them.

Therapy tables Configuration Guide



POSITION 6

Height adjustment and type of control

Chinesport tables can have an electric or hydraulic height adjustment. The option with electric adjustment also has various drive types, depending on the number of electric actuators. This design choice was made for one of the most reliable electric actuators available on the market. The standard mains supply is 220V 50-60Hz and all drives work with very low safety voltage. Clients can make a specific request for 110V 50-60Hz supply and a non-SCHUKO plug. Some models can be powered by a 24V battery that come with a 110/220V 50-60Hz battery charger. The minimum protection rating for all electric components is IP 54.



A Electric, with foot switch

This type of control can be used to activate the electric height adjustment function. The foot switch is placed on the floor. The connection cable allows for complete freedom of movement.



B Electric, with perimeter foot switch

This is an electric height adjustment system that can be reached from anywhere around the table, without having to locate the foot switch. It is activated by simply pressing or lifting the perimeter bar.



C Hydraulic, with double foot switch

This hydraulic adjustment is ideal for situations where there is no power supply or when you need to move the table easily.



D Electric, with hand control
Innovative ergonomic hand control with magnet. Users can attach it wherever they wish on the table frame. It comes with connection cable. The series Unix, Alux and Mini Top have a different hand control.



M Electric with hand control and keypad
This option is exclusive to ALUX tables for beauty and massage. It is a dual electric adjustment for height and Trendelenburg / anti-Trendelenburg position. In particular, this configuration features a hand control plus another control panel on the base frame on both long sides.



E Electric, two motors with foot switch
It can be used with models that have two electric actuators, one for table height adjustment and one for the seat or backrest section.



F Electric, two motors with hand control
This option has an ergonomic hand control with magnet that users can attach wherever they wish on the table frame. This is an alternative to the foot switch control for models that have two electric actuators, one for table height adjustment and one for the seat or backrest section.



G Electric, with dual control
This is an alternative option for electric height adjustments. The control is securely placed in a special cradle, which can be attached at the base of the structure if so required by the user. The control can also be used as hand control, or attached to the table frame by using the magnet.



Therapy tables Configuration Guide



POSITION 7

Other options

Various additional options are available to configure our Chinesport tables after having chosen the height adjustment system and the type of drive. At the moment, we can add wheels to move the table easily, and safety siderails or armrests to provide greater comfort for the patient. These options can also be combined. We would be grateful if you could provide us with any preferred options by writing directly to our International Sales Department at the address export@chinesport.it

- 1 No options**
When inserting the product configuration code, please specify if you are interested in any of the options available for position 7 of the code.



- 2 With wheels for transportation**
Swiveling, non-marking wheels with a central braking system activated via a single foot pedal. Moving the table is extremely comfortable and requires no effort. The table cannot be moved when a patient is lying on the table.




- 3 Safety siderails**
Siderails offer greater safety to patients while they are lying on the table, or when the professional needs to step away temporarily. They can be fitted as and when needed thanks to an easy manual mechanism.



- 4 With wheels and safety siderails**
Please indicate if you would like the table to be fitted with wheels and safety siderails, for greater patient safety. The wheels braking system is controlled by a single foot pedal.



Here is an outline of the height-adjustable tables from the VISIT range for examination and manipulation. The movable sections are adjusted via a system that is servo-assisted by one or more gas springs.




5 Additional armrests
Padded supports are available in the same color as the upholstery of the chosen model. They can be rotated by 180° and moved by 12 cm. They are designed to achieve the patient postures also specified under the Global Postural Re-education method. (1)



6 With wheels and armrests
Here you can select to add wheels to the table configuration, in addition to the armrests that can rotate by 180°.



7 Compact safety siderails
These are alternative safety siderails that can be attached exclusively to the VISIT TRANSFER table. A characteristic of these siderails is that they can be attached or removed without adding any bulk at the sides.



(1) RPG Method - Global Postural Re-education Method

Therapy tables Configuration Guide



POSITION 8

Type of synthetic leather design

- W Seamless rounded edges
- X Squared edges



COMFORT AND FLEXIBILITY



TEAR-RESISTANT



WATERPROOF



EASY MAINTENANCE



FIRE-RESISTANT



RESISTANT TO UV RADIATIONS



FUNGICIDAL TREATMENT

GiottoCollection Standard Range

Giotto is a coated fire-resistant fabric. It is sturdy and easy to care, making it particularly suitable for use in hospitals, care homes, spas, rehabilitation centers and gyms. Giotto has a leather effect surface and shiny finish, and comes in a wide array of colors. Giotto contains Vinyzene, an anti-microbial that is added during the processing stages to create a protective barrier against dust mites, bacteria and fungi, preventing any bad odors and allergies.

TECHNICAL FEATURES	UNITS	STANDARD	VALUES
Weight	g/m ²	-	550
Height	cm	-	140
Composition	-	-	88% PVC - 12% PES
Abrasion resistance	cycles	as per ISO 5470-2	100.000
Color fastness to dry rubbing	grey scale	UNI EN ISO 105-X12	4/5
Color fastness to wet rubbing	grey scale	UNI EN ISO 105-X12	4/5
Color fastness to light	blue scale	UNI EN ISO 105-B02	6
Color fastness to sweat	grey scale	ISO 105-E04	4
Adhesion	N	UNI EN ISO 2411	60
Elongation at break	N	UNI 1421	long. 450 / 190 across
Resistance to tearing	N	UNI 4818-9 met.B	long. 42 / 35 across

FIRE RESISTANCE UNI 9175 class 1 IM - EN 1021-2 BS 5852-1

PADDING DENSITY AND DEPTH


Chinesport normally uses a padding depth and density that is 40mm and 30kg/m³, except for the Bobath range. However, our Upholstery Department is able to deal with particular non-standard requests. In such cases, an additional charge may be applied.

PANEL STABILITY

Panels are secured to the frame by means of bushings. If required, these bushings can also be easily replaced when the padded section wears out or is damaged.

POSITION 9

Color of section cover.



CHOOSING THE UPHOLSTERY COLOR:

Upon placing an order, please remember to specify the upholstery color code. Choose among a standard range of 24 colors.

A WHITE	N IVORY	8 BEIGE	7 ASHGREY	K GREY	S DARK BROWN
B CREAM	4 YELLOW	T GREEN LEMON	1 DARK GREEN	6 LIGHT GREEN	E AQUAMARINE
Z ORANGE	G CORAL	F RED	H BURGUNDY	9 FUCHSIA	Q LILAC
R LIGHT CORNFLOWER BLUE	2 LIGHT BLUE	3 MARINE BLUE	L DARK BLUE	M ANTHRACITE	P BLACK



02023 CLEANING SET
To clean and disinfect the padded synthetic leather surface of any examination and treatment table. The set includes 3 x 500 ml bottles. Another set of protection cream is included, 3 jars.




XCAT0020131 COLOR SAMPLES
Upon request from design companies or individual professionals who would like to check the actual synthetic leather color before making a purchase, we are happy to send our samples.


Accessories Cushions for posture

CUSHIONS FOR POSTURE


Postural cushions are available in various shapes and sizes to meet the professional's varied requirements to position a patient correctly. Our cushions are made from variable density foam and have a synthetic leather cover which is fire-resistant and anti-bacterial, easy to clean and durable. Standard density is 21 kg/m³. Cushion covers have a zip for easy removal, and come in a choice of upholstery from an assortment of 24 color options. Cushions can be purchased in the same type of upholstery as the examination or therapy table.




AC0024.W ? FACE CUSHION
L 31 x D 23 x H 6




AC0025.W ? CUSHION WITH HOLE
L 40 x D 34 x H 9





10321.W ? LORDOSIS SUPPORT
L 29 x D 14 x H 5



10300.W ? SEMI-CYLINDER 1 - L 48 x D 25 x H 18
10310.W ? SEMI-CYLINDER 2 - L 60 x D 40 x H 18



09940.W ? CUBE 1 L 40 x D 40 x H 40
09950.W ? CUBE 2 L 50 x D 50 x H 50





09700.W ? RECTANGLE 1 - L 40 x D 20 x H 10
09730.W ? RECTANGLE 2 - L 40 x D 30 x H 5
09740.W ? RECTANGLE 3 - L 40 x D 30 x H 10
09820.W ? RECTANGLE 4 - L 80 x D 40 x H 40

THERAPY TABLES



11250.W ? LEG CUSHION
L 77 x D 50 x H 20/15



MADE-TO-MEASURE CUSHIONS


On request, we can manufacture cushions with non-standard size and density.

Feel free to send us your request by e-mail to our address: export@chinesport.it

KEY: W = width; D = depth; H = height; all dimensions are in centimeters (cm).

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.






AC0881.W ? FACE CUSHION 2


A cushion with a special shape and density which makes comfortable the prone position of the patient during treatment and therapies. This cushion is advisable for offering great relaxation to patients. Dimensions: W 29 x D 31 x H 9 cm.

REMOVABLE COVERS


All Chinesport mattresses and cushions have a synthetic leather cover which is fire-resistant and anti-bacterial, highly durable, removable and washable.



- 09600.W ? WEDGE 1 - W 25 x D 25 x H 10
- 09610.W ? WEDGE 2 - W 35 x D 35 x H 10
- 09620.W ? WEDGE 3 - W 60 x D 45 x H 15
- 09630.W ? WEDGE 4 - W 60 x D 45 x H 30
- 09640.W ? WEDGE 5 - W 60 x D 60 x H 15
- 09650.W ? WEDGE 6 - W 60 x D 60 x H 20



- 10000.W ? CYLINDER 1 - W 20 x Ø 8
- 10010.W ? CYLINDER 2 - W 35 x Ø 15
- 10020.W ? CYLINDER 3 - W 50 x Ø 15
- 10030.W ? CYLINDER 4 - W 50 x Ø 25
- 10050.W ? CYLINDER 5 - W 100 x Ø 25



- 10130.W ? RIGID ROLL 1 - W 100 x Ø 30
- 10150.W ? RIGID ROLL 2 - W 100 x Ø 40
- 10160.W ? RIGID ROLL 3 - W 100 x Ø 50
- 10200.W ? RIGID ROLL 4 - W 50 x Ø 30



Cushions come in various sizes, they are elastic and mould easily but always return to their original shape. They offer an ideal support to the lumbar and lateral region, as well as to the upper and lower limbs. They are covered in a washable, breathable fabric.


- 01449 POSTURAL CUSHION 100 - ø cm 19 x 17
- 01450 POSTURAL CUSHION 110 - ø cm 19 x 31
- 01451 POSTURAL CUSHION 120 - ø cm 25 x 37
- 01452 POSTURAL CUSHION 8 - cm 40 x 20 x 13 H
- 01453 POSTURAL CUSHION 88 - cm 50 x 28 x 13 H

01426 POSTURAL SET

This set includes 5 cushions, one for each type specified above, and a carry bag to store the items, size 33 x 33 x 56 h cm.

Accessories for height adjustable therapy tables

ELECTRICAL ACCESSORIES



AC0037 ADDITIONAL FOOT PEDAL 1
For models fitted with a single motor. The additional foot pedal provides another way to control the table height adjustment. It can be used on both sides.



AC0038 ADDITIONAL FOOT PEDAL 2
For models fitted with two motors. The additional foot pedal provides another way to control the table height adjustment, as well as the motorized sections. This accessory can be used on both sides.



AC0039 ADDITIONAL HAND CONTROL 1
For models fitted with a single motor. The additional hand control provides another way to control the table height adjustment. This hand control has a magnet so it can be easily attached to the table frame, without any restrictions on the location.




AC0040 ADDITIONAL HAND CONTROL 2
For models fitted with two motors. The additional hand control provides another way to control the table height adjustment, as well as to adjust the motorized sections. This hand control has a magnet so it can be easily attached to the table frame, without any restrictions on the location.



AC0083 ADDITIONAL CONTROL 1
It can be used both as a foot pedal in combination with the corresponding holder, code AC0997, or as a hand control. It has a magnet that enables it to be attached anywhere on the frame.



AC0997 FOOT SWITCH HOLDER
Can be applied on either side of the base frame as preferred by the carer, and enables the accessory AC0083 to be fitted. Not compatible for therapy tables with perimeter footswitch (option B).



AC0578 SAFETY DEVICE
This device disconnects the mains supply from the height adjustment system. This stops any inappropriate use when the carer is away. The device consists of a manual switch.



AC0017 AUXILIARY BATTERY
Enables the table to be used without power supply. It can be attached to the table by plugging the battery cable directly into the actuator. The battery comes with a support so it can be secured to the frame base. Works for 10 height-adjustment cycles. Power supply: 1.2Ah 24V .

FOR PROTECTION AND HYGIENE



AC0031 PAPER ROLL HOLDER F

A metal support with adjustable width from 50 cm to 65 cm that is attached to the upper part of the couch. Recommended for use with paper rolls from set AC0034. This support can be attached to the foot end of models from the Unix, Visit and Ther Series, or to the head end of models from the Test Series.



AC0699 PAPER ROLL HOLDER H

A metal support with adjustable width having similar characteristics to support AC0031, but designed to be attached only to the head end of models from the Unix, Visit, Ther and Test Series. In particular, for models from the Test Series it is compatible for use with push handle AC0698. Available for therapy tables in production from September 2015 only.



AC0032 PAPER ROLL HOLDER G

A metal support that is attached to the foot end of the couch base. It allows for the paper roll to be placed vertically if required. This support can be used with rolls that are up to 68 cm long. This support can be attached to the foot end of models from the Unix, Visit and Ther Series, or to the head end of models from the Sinthesi Series.



AC0034 STANDARD PAPER ROLL SET

The set includes 6 x 60 cm wide rolls. Each roll is 85 meters long.



AC0035 MAXI PAPER ROLL SET

The set includes 2 x 90 cm wide rolls. Each roll is 100 meters long.



AC0036 PROTECTIVE SHEET SET

This set contains 10 synthetic fabric sheets with elasticated border. This sheet can be used on any examination and treatment table.

PADDED SUPPORTS AND ELEMENTS



AC0020.W BREATHING HOLE PLUG

For any examination, massage and treatment tables that provide for the option of a standard nose/mouth breathing hole. When ordering, please specify the chosen upholstery color (can be supplied individually).



AC0307 CYLINDER HOLDER

It can only be attached to model VISIT TRANSFER code LV351. It is fitted to the long side of the table base frame (can be supplied individually).



AC0021 RPG BAR FOR COUCH

This bar is height-adjustable and is fitted with a rope and strap AC0056. This bar can also be used for exercises from the RPG (Global Postural Re-education) method. It can be attached to couches from the THER Series and models UNIX 2, EVO UNIX 2, EVO UNIX PLUS.



AC0018.W LATERAL ARMRESTS

They are used to extend the upper body section for larger individuals, or to improve comfort when positioning the arms. Length is 58 cm and width is 15 cm. (can be supplied in combination).



AC0313 DRIP STAND SUPPORT

To be used with the drip stand accessory code AC0027. It can only be attached to models from the VISIT series (can be supplied in combination).



AC0022 RPG WALL BAR

This bar has the same characteristics and intended use as accessory AC0021, i.e. it is meant for exercises from the RPG (Global Postural Re-education) method. It is attached to the wall and can therefore be used with any examination or therapy table.

Accessories for height adjustable therapy tables

SPECIFIC ACCESSORIES FOR TEST SERIES



AC0026 FASTENING DEVICE

This device can only be applied to models from the Test Series, 63 cm wide. It is used to configure the table/armchair with various types of armrests and other accessories.



AC0027 DRIP STAND

A chromed steel drip stand that is attached to models from the Visit Series. To be ordered together with the pair of supports accessory code AC0313. If attached to models from the Test series, please use the fastening device AC0026.



AC0028 ARMCHAIR ARMREST

This armrest is made from thermoplastic material that can be easily disinfected. When fitted into the fastening device AC0026, it can be easily adjusted in height and rotation. It can be attached in a different position from the other side.



AC0029 PHLEBOTOMY ARMREST

This phlebotomy armrest is made from thermoplastic material that can be easily disinfected. When fitted into the fastening device AC0026, it can be easily adjusted in height and rotation. It can be attached in a different position from the other side in a different position from the other side.




AC0710.W ? PADDED ARMREST

Padded, synthetic leather armrest for greater comfort. When fitted into the fastening device AC0026, it can be easily adjusted in height and rotation. Please select the upholstery color.



AC0698 PUSH HANDLE

It can only be applied to model Test Mobil code LE 211 to facilitate transportation when the patient is lying on the table.



AC0023.W ? CERVICAL CUSHION

W 30 x D 15 x H 8,5



AC0545.W ? PLUG FOR GYNECOLOGICAL EXAMINATION TABLE

This plug allows the table to be used for gynecological or urological purposes - ref. models LE151 or LE 351 - and the patient can either lie supine or prone. You can select a color that matches the couch upholstery, so please remember to specify the upholstery color along with the accessory code.




? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.



Visit Line



This is the classic range of examination and treatment tables consisting mainly of two sections, with a 68 cm-long head section that can be used as a backrest. The backrest can have positive and negative tilt adjustment: up to +70° /- 20°. Height is adjusted via one of the most reliable motors available on the market. A hydraulic version is also available as an alternative. The frame is unique for its simplicity and harmonious design, which is protected by registered Community design. Some of our tables are also available with larger wheels so they can be used to transport patients. This configuration can be completed by fitting a set of safety siderails.




"I bought a VISIT 1 couch four years ago and I use it several times a week. I had seen Chinesport couches in several hospital and really liked them. I am absolutely pleased not only with its features, but also because of the numerous positive comments that I receive from my patients, particularly the elderly ones, and for its ease of access. It allows me to do away with those inconvenient and dangerous steps." Prof. Carlo Moreschi - Legal Medical Examiner (Udine - Italy)




Pos. 3 - 4 - 5 Type of sections


111 VISIT 1



121 VISIT LARGE



131 VISIT EXTRA LARGE



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 3** Safety siderails
- 4** Wheels and siderails


Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Model shown - code LV 111 A1 W6



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+70° -20°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	70 / 74

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0020.W? BREATHING HOLE PLUG
- AC0881.W? FACE CUSHION 2
- AC0313 DRIP STAND SUPPORTS **
- AC0027 DRIP STAND **

CUSHIONS FOR POSTURE:

Refer to the first pages

* Foot switch holder not compatible if perimeter footswitch (option B) is chosen


** Accessories AC0313 and AC0027 are not available for model LV131.

upholstery
color code



Pos. 3 - 4 - 5 Type of sections

141 VISIT NAR



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 3** Safety siderails
- 4** Wheels and siderails

Pos. 8 Synthetic leather design

- W** Seamless rounded edges


Pos. 9 Section upholstery color

24 options available



Detail of electrical height adjustment with perimeter foot switch - option B

For further information on options and accessories, please refer to the first pages.



Model shown - code LV 141 A1 W6

Configuration example for a treatment table from the VISIT range.



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+70° -20°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	70

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0020.W? BREATHING HOLE PLUG
- AC0313 DRIP STAND SUPPORTS
- AC0027 DRIP STAND

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code



Pos. 3 - 4 - 5 Type of sections

151 VISIT TREND



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges


Pos. 9 Section upholstery color

24 options available



Detail of hydraulic height adjustment with double foot switch - option C.

For further information on options and accessories, please refer to the first pages.



Model shown - code LV 151 A1 W6



Configuration example for a treatment table from the VISIT range.



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+70° -20°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	70

ACCESSORIES:

ELECTRICAL ACCESSORIES:
 AC0037 ADDITIONAL FOOT PEDAL 1
 AC0039 ADDITIONAL HAND CONTROL 1
 AC0083 ADDITIONAL CONTROL 1
 AC0578 SAFETY DEVICE
 AC0017 AUXILIARY BATTERY
 AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:
 AC0031 PAPER ROLL HOLDER F
 AC0699 PAPER ROLL HOLDER H
 AC0032 PAPER ROLL HOLDER G
 AC0034 STANDARD PAPER ROLL SET
 AC0035 MAXI PAPER ROLL SET
 AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:
 AC0020.W? BREATHING HOLE PLUG
 AC0313 DRIP STAND SUPPORTS
 AC0027 DRIP STAND

CUSHIONS FOR POSTURE:
 Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code



Pos. 3 - 4 - 5 Type of sections

161 VISIT FLEXION



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 3** Safety siderails
- 4** Wheels and siderails

Pos. 8 Synthetic leather design

- W** Seamless rounded edges


Pos. 9 Section upholstery color

24 options available



Detail of electrical height adjustment with dual control - option G.

For further information on options and accessories, please refer to the first pages.



Model shown - code LV 161 A1 W6



TECHNICAL FEATURES

Number of sections	4
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+70° -20°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	70

ACCESSORIES:

ELECTRICAL ACCESSORIES:
 AC0037 ADDITIONAL FOOT PEDAL 1
 AC0039 ADDITIONAL HAND CONTROL 1
 AC0083 ADDITIONAL CONTROL 1
 AC0578 SAFETY DEVICE
 AC0017 AUXILIARY BATTERY
 AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

AC0031 PAPER ROLL HOLDER F
 AC0699 PAPER ROLL HOLDER H
 AC0032 PAPER ROLL HOLDER G
 AC0034 STANDARD PAPER ROLL SET
 AC0035 MAXI PAPER ROLL SET
 AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

AC0020.W? BREATHING HOLE PLUG
 AC0313 DRIP STAND SUPPORTS
 AC0027 DRIP STAND

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code




Pos. 3 - 4 - 5 Type of sections


211 VISIT MOBIL




212 VISIT MOBIL BASIC



221 VISIT MOBIL LARGE



222 VISIT MOBIL LARGE BASIC



Pos. 6 Height adjustment

- A** Electric, with foot switch (*only LV211, LV221*)
- C** Hydraulic, with double pedal
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 3** Safety siderails


Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Model shown - code LV 211 D1 W6



Treatment table configuration with safety siderails - option 3.

TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+70° -20°
Wheels diameter (mm)	120
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	70 / 74

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0020.W? BREATHING HOLE PLUG
- AC0313 DRIP STAND SUPPORTS
- AC0027 DRIP STAND

CUSHIONS FOR POSTURE:


Refer to the first pages

upholstery
color code





Pos. 3 - 4 - 5 Type of sections

351 VISIT TRANSFER



352 VISIT TRANSFER BASIC



Model shown - code LV 351 C1 W6

Pos. 6 Height adjustment

- C Hydraulic, with double foot switch
- D Electric, with hand control
- G Electric, with dual control

Pos. 7 Other options

- 1 No options
- 7 Compact safety siderails

Pos. 8 Synthetic leather design

- W Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



Compact safety siderails - option 7.

For further information on options and accessories, please refer to the specific pages.

Configuration example for a treatment table from the VISIT range.



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	200
Height adjustment	Electric (battery) / Hydraulic
Minimum/maximum table height (cm)	55-97
Mains supply	N/A - Battery charger 110-220V, 50-60 Hz
Actuator thrust (N)	6000
IP Protection	X4
Head section adjustment	Gas spring-assisted
Head section tilt	+70° -20°
Wheels diameter (mm)	150
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	86

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0071 SPARE BATTERY IPX5
- AC0318 WALL BATTERY CHARGER
- AC0997 FOOT SWITCH HOLDER
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0034 STANDARD PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET
- AC0307 CYLINDER HOLDER

PADDED SUPPORTS AND ELEMENTS:

- AC0020.W? BREATHING HOLE PLUG
- AC0313 DRIP STAND SUPPORTS
- AC0027 DRIP STAND

CUSHIONS FOR POSTURE:

Refer to the first pages

upholstery color code



The quality of Chinesport products is the result of ongoing efforts for improvement and highest attention to detail. For example, the central braking system with one single, ergonomic foot switch control can be activated easily and noiselessly.




Here is an outline of the height-adjustable tables from the VISIT range for examination and manipulation. The movable sections are adjusted via a system that is servo-assisted by one or more gas springs. The frame is sturdy and stable, with a design that is protected by registered Community design.

TherLine




A major difference from the VISIT range is given by the head section length which is only 52 cm. For this reason, it cannot be used as a backrest, but it can be fitted with armrests to provide greater comfort for the patient when in prone position. These tables were developed for physiotherapy, osteopathy or for the so-called "Global Postural Re-education" method, and they can be configured with up to six sections to achieve different postural positions for the patient. There is also a model that can be arranged in a seated position for particular use in gynecology. Another difference from the VISIT range is that the head section has a greater tilt adjustment, i.e. positive tilt up to $+85^\circ$ and negative tilt up to -35° , plus the nose/mouth breathing hole is standard.




Pos. 3 - 4 - 5 Type of sections

111 THER 1



121 THER PLUS 1



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double pedal
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.



Model shown - Code LH 111 B2 WR



Folding lateral armrest, accessory AC0018.W

TECHNICAL FEATURES

Number of sections	2 / 4
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	76 / 78

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0997 FOOT SWITCH HOLDER *
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0018.W? LATERAL ARMRESTS
- AC0021 RPG BAR FOR COUCH

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code




Pos. 3 - 4 - 5 Type of sections

112 THER 2



122 THER PLUS 2



Model shown - code LH 122 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double pedal
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges


Pos. 9 Section upholstery color

24 options available

Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.



Example of portioned trendelenburg position.

TECHNICAL FEATURES

Number of sections	3 / 5
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	80 / 82

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0997 FOOT SWITCH HOLDER *
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

ADDITIONAL SUPPORTS:

- AC0018.W? LATERAL ARMRESTS
- AC0021 RPG BAR FOR COUCH

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code





Pos. 3 - 4 - 5 Type of sections

113 THER NAR



123 THER PLUS NAR



Model shown - CODE LH 123 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



TECHNICAL FEATURES

Number of sections	3 / 5
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	80 / 82

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0997 FOOT SWITCH HOLDER *
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0018.W? LATERAL ARMRESTS
- AC0021 RPG BAR FOR COUCH

CUSHIONS FOR POSTURE:

Refer to the first pages

Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code





Pos. 3 - 4 - 5 Type of sections

114 THER FLEXION



124 THER PLUS FLEXION



Model shown - code LH 124 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



TECHNICAL FEATURES

Number of sections	4 / 6
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	82 / 84

ACCESSORIES:

- ELECTRICAL ACCESSORIES:**
- AC0037 ADDITIONAL FOOT PEDAL 1
 - AC0039 ADDITIONAL HAND CONTROL 1
 - AC0997 FOOT SWITCH HOLDER *
 - AC0083 ADDITIONAL CONTROL 1
 - AC0578 SAFETY DEVICE
 - AC0017 AUXILIARY BATTERY
- FOR PROTECTION AND HYGIENE:**
- AC0031 PAPER ROLL HOLDER F
 - AC0699 PAPER ROLL HOLDER H
 - AC0032 PAPER ROLL HOLDER G
 - AC0034 STANDARD PAPER ROLL SET
 - AC0035 MAXI PAPER ROLL SET
 - AC0036 PROTECTIVE SHEET SET

- PADDED SUPPORTS AND ELEMENTS:**
- AC0018.W? LATERAL ARMRESTS
 - AC0021 RPG BAR FOR COUCH

- CUSHIONS FOR POSTURE:**
Refer to the first pages

* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

Configuration example for a treatment table from the VISIT range.





For further information on options and accessories, please refer to the first pages.

upholstery color code



Pos. 3 - 4 - 5 Type of sections

118 THER GYN



Model shown - CODE LH 118 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.

Padded thigh holders are included.



Sitting position with accessory headrest code AC0023.W

TECHNICAL FEATURES

Number of sections	3
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	80

ACCESSORIES:

ELECTRICAL ACCESSORIES:
 AC0037 ADDITIONAL FOOT PEDAL 1
 AC0039 ADDITIONAL HAND CONTROL 1
 AC0997 FOOT SWITCH HOLDER *
 AC0083 ADDITIONAL CONTROL 1
 AC0578 SAFETY DEVICE
 AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:
 AC0031 PAPER ROLL HOLDER F
 AC0032 PAPER ROLL HOLDER G
 AC0034 STANDARD PAPER ROLL SET
 AC0035 MAXI PAPER ROLL SET
 AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:
 AC0021 RPG BAR FOR COUCH

CUSHIONS FOR POSTURE:
 Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code





Pos. 3 - 4 - 5 Type of sections

117 THER MODULAR



127 THER PLUS MODULAR



Model shown - CODE LH 127 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.

Configuration example for a treatment table from the VISIT range.



TECHNICAL FEATURES

Number of sections	3 / 5
Max safety load (kg)	200
height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	82 / 84

ACCESSORIES:

- ELECTRICAL ACCESSORIES:**
 AC0037 ADDITIONAL FOOT PEDAL 1
 AC0039 ADDITIONAL HAND CONTROL 1
 AC0997 FOOT SWITCH HOLDER *
 AC0083 ADDITIONAL CONTROL 1
 AC0578 SAFETY DEVICE
 AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
 AC0699 PAPER ROLL HOLDER H
 AC0032 PAPER ROLL HOLDER G
 AC0034 STANDARD PAPER ROLL SET
 AC0035 MAXI PAPER ROLL SET
 AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0018.W? LATERAL ARMRESTS
 AC0021 RPG BAR FOR COUCH

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code




Pos. 3 - 4 - 5 Type of sections

109 THER OSTEO D



135 THER OSTEO



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.



Model shown - CODE LH 135 B2 WZ



AC0019.W? **LOAD-RELIEF ARMREST**

The head section can be turned by up to 180°.

TECHNICAL FEATURES

Number of sections	1 / 2
Max safety load (kg)	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted (only LH135)
Head section tilt	+85° -35°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	50
Weight (kg)	70

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0997 FOOT SWITCH HOLDER *
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H (only LH135)
- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0019.W? LOAD-RELIEF ARMREST
- AC0021 RPG BAR FOR COUCH

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code



Pos. 3 - 4 - 5 Type of sections

126 THER PLUS GLOBAL



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control



Model shown - CODE LH 126 B2 WZ



Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



The two lateral armrests can be turned by up to 180° and moved by 12 cm.

TECHNICAL FEATURES

Number of sections	6
Max safety load (kg)	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	50
Weight (kg)	76

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0997 FOOT SWITCH HOLDER *
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H
- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0021 RPG BAR FOR COUCH
- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:

Refer to the first pages

Configuration example for a treatment table from the VISIT range.




For further information on options and accessories, please refer to the first pages.

upholstery color code



Global Postural Re-education, more frequently identified with the acronym RPG, originates from the book “Le Champ Clos” (The Closed Field) written by Professor Ph. E. Souchard in 1981. The end goal of RPG is to globally restore the body’s static and dynamic balance after it had been compromised due to spontaneous defense mechanisms that were kick-started by external aggressions. The principle of RPG is based on postures involving an active muscular stretching. These postures tend to highlight the muscle tension chains in tendons and membranes around muscles (which lock the individual in an unbalanced position). These postures create an active elongation (active participation from the individual) and a subsequent relaxation of the muscle chains that were affected by spasm and stiffness. Unlike other types of treatment, these postures need a systematic, global approach in order to be effective. Special focus needs to be placed on breathing. There are countless benefits that can be felt in many areas: motor, static, digestive, circulatory, mental and emotional. An overall balance requires the respiratory muscles, and particularly the diaphragm, to work efficiently.

The therapy table THER PLUS GLOBAL allows us to correctly perform the various postures used in the Global Postural Re-education method. In particular, the two lateral armrests can be turned by up to 180° and moved by 12 cm.



uniX Line




This range is an alternative to the classic ranges called VISIT and THER used for examinations, treatments and physiotherapy. In particular, the distinguishing feature in these tables is that height adjustment occurs while the individual sections remaining completely motionless. This makes them particularly suitable for smaller environments or to make the best use of available space, as well as to maintain the position chosen for exercises once the pulley therapy system has been set up. You can opt for a simplified base frame and configure the head section with a backrest that is 68 cm long, instead of a smaller, 52 cm long section that can come with armrests. The second type of head section can have positive and negative tilt adjustment: up to $+85^\circ$ / -35° .




Pos. 3 - 4 - 5 Type of sections

111 UNIX 1



112 UNIX 2



Pos. 6 Height adjustment

D Electric, with hand control

Pos. 7 Other options


- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

W Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



Model shown - code LC 111 D2 W7



Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.

TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	170
Height adjustment	Electric
Minimum/maximum table height (cm)	55-93
Mains supply	230 V 50/60 Hz 1,5A
Actuator thrust (N)	6000
IP Protection	44
Head section adjustment	Gas spring-assisted
Head section tilt - Unix 1	+70°
Head section tilt - Unix 2	+85° -35°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	70

ACCESSORIES

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET
- AC0699 PAPER ROLL HOLDER H

PADDED SUPPORTS AND ELEMENTS:

- AC0020.W? BREATHING HOLE PLUG (Unix 1 only)
- Unix 2 comes with breathing hole plug as standard
- AC0021 RPG BAR FOR COUCH (Unix 2 only)
- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:


Refer to the first pages

upholstery color code



Pos. 3 - 4 - 5 Type of sections

121 EVO UNIX 1



Model shown - code LC 121 B2 W7

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg) - Electric	250
Max safety load (kg) - Hydraulic	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	50-97
Mains supply	230 V 50/60 Hz 2A
Actuator thrust (N)	8000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+70°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	85

ACCESSORIES

FOR PROTECTION AND HYGIENE:
 AC0031 PAPER ROLL HOLDER F
 AC0034 STANDARD PAPER ROLL SET
 AC0035 MAXI PAPER ROLL SET
 AC0036 PROTECTIVE SHEET SET
 AC0699 PAPER ROLL HOLDER H

PADDED SUPPORTS AND ELEMENTS:
 AC0020.W? BREATHING HOLE PLUG
 AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:
 Refer to the first pages

Configuration example for a treatment table from the VISIT range.




For further information on options and accessories, please refer to the first pages.

upholstery color code




Pos. 3 - 4 - 5 Type of sections


122 EVO UNIX 2




123 EVO UNIX PLUS



124 EVO UNIX 2 TREND



125 EVO UNIX PLUS TREND



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- C** Hydraulic, with double foot switch

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Model shown - code LC123B2W7



TECHNICAL FEATURES

Number of sections	2 / 4
Max safety load (kg) - Electric	250
Max safety load (kg) - Hydraulic	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	50 - 97
Mains supply	230 V 50/60 Hz 2A
Actuator thrust (N)	8000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+85° -35°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	85 (LC122; LC123) 90 (LC124; LC125)

ACCESSORIES

FOR PROTECTION AND HYGIENE:
AC0031 PAPER ROLL HOLDER F
AC0034 STANDARD PAPER ROLL SET
AC0035 MAXI PAPER ROLL SET
AC0036 PROTECTIVE SHEET SET
AC0699 PAPER ROLL HOLDER H

PADDED SUPPORTS AND ELEMENTS:
AC0021 RPG BAR FOR COUCH
AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:
 Refer to the first pages

upholstery
color code




Sinthesi

Our Sinthesi couches allow therapists to effortlessly place patients in countless myofascial stretching postures that can be easily maintained for longer periods of time and without strain. In this relaxed position, various techniques can be applied to facilitate recovery in myofascial elongation, increase joint movement range and achieve an antalgic posture as a relief from acute pain. Also, Sinthesi tables provide the ideal starting position for exercises to stimulate central activation and for isometric muscular work.





Pos. 3 - 4 - 5 Type of sections

371 SINTHESI SINGLE



381 SINTHESI PLUS SINGLE



Model shown - code LS 371 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 5** Additional lateral armrests
- 6** Wheels and lateral armrests

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Additional control panel for electric adjustments as standard.

TECHNICAL FEATURES

Number of sections	3 / 5
Max safety load (kg)	180
Height adjustment	Electric
Minimum/maximum table height (cm)	48-100
Mains supply	220V 50/60 Hz - 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Backrest adjustment	Electric
Head section tilt	+54° -52°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	95

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0784 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code




Pos. 3 - 4 - 5 Type of sections

375 SINTHESI NAR



385 SINTHESI PLUS NAR



Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 5** Additional lateral armrests
- 6** Wheels and lateral armrests


Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Model shown - code LS 375 B2 WR



Additional control panel for electric adjustments as standard.

TECHNICAL FEATURES

Number of sections	3 / 5
Max safety load (kg)	180
Height adjustment	Electric
Minimum/maximum table height (cm)	48-100
Mains supply	220V 50/60Hz - 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Backrest adjustment	Electric
Head section tilt	+54° -52°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	95

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0784 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code





Pos. 3 - 4 - 5 Type of sections

372 SINTHESI DOUBLE



382 SINTHESI PLUS DOUBLE



Model shown - code LS 372 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 5** Additional lateral armrests
- 6** Wheels and lateral armrests

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Additional control panel for electric adjustments as standard.

TECHNICAL FEATURES

Number of sections	4 / 6
Max safety load (kg)	180
Height adjustment	Electric
Minimum/maximum table height (cm)	48-100
Mains supply	220V 50/60 Hz - 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Backrest adjustment	Electric
Head section tilt	+54° -52°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	97

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0784 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code





Pos. 3 - 4 - 5 Type of sections

373 SINTHESI SPLIT



383 SINTHESI PLUS SPLIT



Model shown - code LS 383 B2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 5** Additional lateral armrests
- 6** Wheels and lateral armrests

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Additional control panel for electric adjustments as standard.

TECHNICAL FEATURES

Number of sections	4 / 6
Max safety load (kg)	180
Height adjustment	Electric
Minimum/maximum table height (cm)	48-100
Mains supply	220V 50/60Hz - 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Backrest adjustment	Electric
Head section tilt	+54° -52°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	97

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0784 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:

Refer to the first pages


* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code





Pos. 3 - 4 - 5 Type of sections

374 SINTHESI MITO



384 SINTHESI PLUS MITO



Model shown - code LS 374 A2 WR

Pos. 6 Height adjustment

- A** Electric, with foot switch
- B** Electric, with perimeter foot switch
- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control
- G** Electric, with dual control



Additional control panel for electric adjustments as standard

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation
- 5** Additional lateral armrests
- 6** Wheels and lateral armrests

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.

TECHNICAL FEATURES

Number of sections	6 / 8
Max safety load (kg)	180
Height adjustment	Electric
Minimum/maximum table height (cm)	48-100
Mains supply	220V 50/60 Hz - 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Backrest adjustment	Electric
Head section tilt	+54° -52°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	101

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0083 ADDITIONAL CONTROL 1
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY
- AC0997 FOOT SWITCH HOLDER *

FOR PROTECTION AND HYGIENE:

- AC0032 PAPER ROLL HOLDER G
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0784 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:

- AC0022 RPG WALL BAR

CUSHIONS FOR POSTURE:

Refer to the first pages

* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code



Foundation course on myofascial release techniques using Sinthesi therapy tables and global postural analysis.


Learning objectives:

Participants shall learn how myofascial structures adapt to poor postures in order to economize on movements and avoid pain in the joint movement range. A global postural analysis through GPS enables therapists to identify which groups are requiring treatment in order to regain appropriate posture and movement.


This course shall demonstrate how Sinthesi treatment tables have become an integral part of postural training and treatment for many therapists. With 8 independently mobile sections and uninterrupted, noiseless functions of movement control, therapists can offer a new and unique range of myofascial treatments that provide an efficient and effective treatment for many disorders.

Our Sinthesi tables allow therapists to effortlessly place patients in countless myofascial stretching postures that can be easily maintained for longer periods of time and without strain. In this relaxed position, various techniques can be applied to facilitate recovery in myofascial elongation, increase joint movement range and achieve an antalgic posture as a relief from acute pain. Also, Sinthesi tables provide the ideal starting position for exercises to stimulate central activation and for isometric muscular work.

The myofascial release techniques course is now available worldwide from Chinesport. It is the ideal starting point for trainee practitioners striving to improve their skills. Participants will learn how to use the unique MITO - Sinthesi therapy table within a global approach towards myofascial therapy, with the aim to restore the natural myofascial length and improve posture and movement.



This video shows our postural treatment table MiTo PLUS, with an introduction to how it can be adjusted and the positions that can be achieved. (excerpt from Dr. Francesco Mignani's introductory DVD course)



01778


MANUAL FOR MI.TO. TREATMENT TABLE

This document was developed by its author with the aim to offer useful therapeutic indications through the use of a special therapy table for active and passive treatments via postural exercises. Special emphasis is given to respiratory training, awareness of the vertebral column function, mobilization exercises, muscle stretching, decompression of discs and tractions, strengthening and stabilization of the rachis, manual treatment for conditions (e.g. spondylolisthesis, lumbosciatica due to hernia, painful scoliosis, hyperkyphosis) and application in sporting activities. Author Dr. Francesco Mignani; in Italian/English; Size 15 x 21; Pages 30

“Interpreting and understanding any dysfunctional mechanisms caused by a change in the balance of the body unit leads us to consider how important it is to regain such balance, by working on restoring harmony in muscle chains, on the principles of their functional system and their potential correlations (ascending, descending or mixed).”

“If one or more elements, whether internal or external, interfere with this harmony, our body tries to maintain balance by putting the lack of pain first. In order not to suffer pain, our body bends, extends, tilts, reduces its mobility and tries to seek comfort through these defensive adaptations that are certainly less economical. Unfortunately, this leads to significant energy loss particularly for myofascial structures.”
(excerpt from manual)

FEEDBACK ON MI.TO. THERAPY TABLES



Prof. **Benedetto Toso**, Posturology Lecturer at the Faculty of Motor Sciences of the Università Cattolica del Sacro Cuore in Milan, Founder and President of the Back School Association.

"The MITO therapy table is a useful tool for professionals so muscle stretching exercises can be performed effectively and by each muscle group, avoiding any potential compensation. Also, in the event of acute lumbago, it enables the patient to passively get into an antalgic posture that helps to re-balance the rachis group and reduce pain. All these features make it a precious tool to have both in physical therapy practices and gyms. The only limitation (which is not in fact a limitation but an advantage) is that you will need to have more than one therapy table in your gym as it helps patients to perform numerous exercises in the most appropriate way".



Dr. **Francesco Mignani**, with a degree in Physiotherapy, specialized in Osteopathy, Motor Sciences teacher, massage therapist and physiotherapist, Director of the Centro Fisioproject in Santa Maria degli Angeli (Assisi), Italy

- *Our MiTo postural table is a multi-functional tool both for direct manual treatments and within a wider postural rehabilitation approach.*
- *Manipulation and posture, posture and manipulation - this interaction is now possible with our MiTo postural table.*
- *Our MiTo postural table facilitates the therapist's work by giving continuity to manual treatments, as it also provides great help to perform postural exercises.*
- *When there is a disorder, the body's overall organization needs a combination of strategies that our MiTo postural table is able to address with its multiple functions.*
- *Caring for the individual as a whole and addressing any complex issues now has an additional tool: the MiTo postural table.*



01945.DVD

POSTURAL EXERCISES WITH THE MI.TO. THERAPY TABLE.

This video offers a theory and practice training course that offers an overview of the possible manual treatments combined with a wider approach involving postural exercises. Our Mi.To. therapy table can also be used to perform passive treatments on the patient, as well as for passive/active and active work. The aim is to regain flexibility and combat stiffness and aging. The course covers all the topics already dealt with in the study manual - code 01778. The Mi.To. therapy table is part of the national education syllabus on the "Mignani Postural Method" and of the Milanino Back School Association - Italy.

In Italian and English; Duration: over 90 minutes



Bobath Line




Large physical therapy tables that are up to two meters long and have a surface area measuring four square meters. The size and stability characteristics mean they are suitable for neurological treatments, as patients can be guided through a wide range of exercises and sensory experiences, as well as postural adjustments. This range also includes the models used in the unique R.I.C. - Continued Intensive Rehabilitation method. Finally, a special version also allows for use with static supine exercises without having to move a patient to a rehabilitation gym.




BobathLine Therapy tables

Pos. 3 - 4 - 5 Type of sections


311 BOBATH ECO D100



411 BOBATH SD100



412 BOBATH SD120



Pos. 6 Height adjustment

For all models

- A** Electric, with foot switch
- D** Electric, with hand control
- G** Electric, with dual control

For model LB412 only

- C** Hydraulic

Pos. 7 Other options


- 1** No options
- 2** With wheels for transportation
- 3** Safety siderails
- 4** Wheels and siderails

Pos. 8 Synthetic leather design


- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



Model shown - code LB 412 A2 W3



Accessible up to cm 18

TECHNICAL FEATURES

Number of sections	1
Max safety load (kg) - LB311	150
Max safety load (kg) - LB411-LB412	Electric 250 / Hydraulic 150
Height adjustment	Electric / Hydraulic (LB412 only)
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N) - LB311	6000
Actuator thrust (N) - LB411-LB412	10000
IP Protection	54
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	30
Density (kg/m ³)	30
Weight (kg)	90 / 93 / 98



ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE

FOR PROTECTION AND HYGIENE:

- AC0035 MAXI PAPER ROLL SET

CUSHIONS FOR POSTURE:


Refer to the first pages

upholstery
color code





Pos. 3 - 4 - 5 Type of sections

413 BOBATH SD150



414 BOBATH SD200



Model shown - code LB 414 A2 W3

Pos. 6 Height adjustment

- A** Electric, with foot switch
- D** Electric, with hand control
- G** Electric, with dual control

Pos. 7 Other options

- 2** With wheels for transportation
- 4** Wheels and siderails

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	250
Height adjustment	Electric
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	10000
IP Protection	54
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	30
Density (kg/m ³)	30
Weight (kg)	140 / 156



ACCESSORIES:

- ELECTRICAL ACCESSORIES:**
 AC0037 ADDITIONAL FOOT PEDAL 1
 AC0039 ADDITIONAL HAND CONTROL 1
 AC0083 ADDITIONAL CONTROL 1
 AC0578 SAFETY DEVICE

FOR PROTECTION AND HYGIENE:
 AC0035 MAXI PAPER ROLL SET

CUSHIONS FOR POSTURE:
 Refer to the first pages

Configuration example for a treatment table from the VISIT range.




For further information on options and accessories, please refer to the first pages.

upholstery color code




Pos. 3 - 4 - 5 Type of sections


321 BOBATH ECO T100



421 BOBATH S T100



422 BOBATH S T120



Pos. 6 Height adjustment

For all models

- A** Electric, with foot switch
- D** Electric, with hand control
- G** Electric, with dual control

For model LB422

- C** Hydraulic

Pos. 7 Other options


- 1** No options
- 2** With wheels for transportation
- 3** Safety siderails
- 4** Wheels and siderails

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



Model shown - code LB 422 A1 W3



Rapid and silent electric adjustment cm 47-97

TECHNICAL FEATURES

Number of sections	2
Max safety load (kg) - LB321	150
Max safety load (kg) - LB421-LB422	Electric 250 / Hydraulic 150
Height adjustment	Electric / Hydraulic (LB422 only)
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N) - LB321	6000
Actuator thrust (N) - LB421-LB422	10000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	75°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	30
Density (kg/m³)	30
Weight (kg)	90 / 93 / 103



ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE

FOR PROTECTION AND HYGIENE:

- AC0035 MAXI PAPER ROLL SET

CUSHIONS FOR POSTURE:


Refer to the first pages

upholstery color code




Pos. 3 - 4 - 5 Type of sections

432 BOBATH S RIC 120



472 BOBATH S RIC TREND 120



Pos. 6 Height adjustment

For model LB432 only

D Electric, with hand control

For model LB472 only

F Electric, two motors with hand control

Pos. 7 Other options

1 No options

2 With wheels for transportation

3 Safety siderails

4 Wheels and siderails

Pos. 8 Synthetic leather design

W Seamless rounded edges


Pos. 9 Section upholstery color

24 options available

Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages and the Pulley Therapy chapter.



Model shown - code LB 432 A2 W3

Suitable for the



TECHNICAL FEATURES

Number of sections	4
Max safety load (kg)	250
Height adjustment	Electric
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	10000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	90°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	30
Density (kg/m³)	30
Weight (kg)	98 / 108



ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0038 ADDITIONAL FOOT PEDAL 2 *
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0040 ADDITIONAL HAND CONTROL 2 *
- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE

FOR PROTECTION AND HYGIENE:

- AC0035 MAXI PAPER ROLL SET

CUSHIONS FOR POSTURE:

Refer to the first pages

* IMPORTANT NOTE: accessories AC0038 and AC0040 are only available for model LB 472


upholstery color code




Bobath Line Therapy tables

Pos. 3 - 4 - 5 Type of sections


341 BOBATH TN100



351 BOBATH FLEX 100



361 BOBATH TREND 100



Pos. 6 Height adjustment

For model LB341

- A** Electric, with foot switch
- D** Electric, with hand control
- G** Electric, with dual control

For models LB351, LB361

- F** Electric, two motors with hand control

Pos. 7 Other options


- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges


Pos. 9 Section upholstery color

24 options available



Model shown - code LB 351 F2 W3

Configuration example for a treatment table from the VISIT range.



TECHNICAL FEATURES

Number of sections	2 / 4
Max safety load (kg)	150
Height adjustment	Electric
Minimum/maximum table height (cm)	47-97
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	-20° +75°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	30
Density (kg/m)	30
Weight (kg)	85 / 93 / 93

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0997 FOOT SWITCH HOLDER *
- AC0578 SAFETY DEVICE

FOR PROTECTION AND HYGIENE:

- AC0035 MAXI PAPER ROLL SET

CUSHIONS FOR POSTURE:

Refer to the first pages



* Foot switch holder not compatible if perimeter footswitch (option B) is chosen

upholstery color code



Pos. 3 - 4 - 5 Type of sections

492 BOBATH S TILT TABLE 2ET



L	B								
1	2	3	4	5	6	7	8	9	

Model shown - code LB 492 F2 X3

Pos. 6 Height adjustment

F Electric, two motors with hand control

Pos. 7 Other options

2 With wheels for transportation

Pos. 8 Synthetic leather design

X Squared edges

Pos. 9 Section upholstery color

24 options available



Removable footplate assembly

For further information on options and accessories, please refer to the dedicated "Tilt Tables" chapter.



Configuration example for a treatment table from the VISIT range.



TECHNICAL FEATURES

Number of sections	6
Max safety load (kg)	250 hoisting / 190 tilt table
Height adjustment	Electric
Minimum/maximum table height (cm)	55-105
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	10000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	-20° +75°
Wheels diameter (mm)	65
braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	170

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE

FOR PROTECTION AND HYGIENE:

- AC0035 MAXI PAPER ROLL SET

PADDED SUPPORTS AND ELEMENTS:

- AC0002 INTRA-EXTRA SANDALS
- AC0003 LATERAL TRUNK SUPPORTS
- AC0004 HAND GRIPS
- AC0005 FASTENING STRAP
- AC0006 STRAP CUSHION
- AC0308 BOBATH TILT TABLE
- AC0311 LEG SEPARATOR
- AC0312 STRAP GUIDE

CUSHIONS FOR POSTURE:

Refer to the first pages

upholstery color code



Bobath Line for supine positions and treatments



Fitting the footplate assembly for use in supine standers

1) Note: A special safety device is available for a different height or tilt adjustment.

2) Note: An additional safety system is available to lower the table manually in the event of power failure.


This model was designed to help professionals and patients in gym rehabilitation programs, avoiding any waste of time and discomfort during transfers from one type of treatment table to the other.

EXAMPLES OF ACCESSORY APPLICATIONS


THERAPY TABLES




Foot restraining sandals accessory AC0002 (2 pieces)




Strap guide accessory AC0312 (single piece)



Leg separator accessory AC0311 (single piece)



Lateral trunk supports accessory AC0003 (2 pieces)



Hand grips accessory AC0004 (2 pieces)



Work table for static exercises in supine position accessory AC0308 (single piece)

Test Line



This is a range of multi-functional tables for examination and for specific treatments. In particular, they can be used in podiatry, gynecology, urology, cardiology, geriatrics, beauty, phlebotomy or even to transport patients within care facilities. The table height can be adjusted electrically or via hydraulic piston, while the moving sections are servo-assisted by gas spring. In some versions, the seating section is tilt-adjustable via another motor. As a result, these tables can easily change into a seated, relaxation, Trendelenburg and anti-Trendelenburg position even with the patient lying on the table. A series of accessories is also available for specified uses.



For example, this picture shows one configuration that includes the following accessories:

AC0026 FASTENING DEVICE N. 2 PIECES

AC0028 ARMCHAIR ARMREST N. 2 PIECES

Test Line Multi-functional tables

For all Test Series examination, treatment or postural tables, the possible configurations can be integrated by special accessories that are designed for particular uses. The key component for attaching these accessories is a chromed steel fastening device that can be fitted in various positions on both sides of the seat section. Each fastening device code AC0026 allows for attaching one single accessory at a time, which can then be adjusted in height, rotation and tilt.

AC0028 ARMCHAIR ARMREST
Made from polyurethane, can be supplied individually.

EXAMPLE 1
with the addition of:
• n. 2 fastening devices
• n. 2 armchair armrests

AC0710.W? PADDED ARMREST
Made from polyurethane, can be supplied individually.

EXAMPLE 2
with the addition of:
• n. 1 fastening device
• n. 1 padded armrest

AC0029 PHLEBOTOMY ARMREST
Made from polyurethane, can be supplied individually.

EXAMPLE 3
with the addition of:
• n. 2 fastening devices
• n. 1 armchair armrest
• n. 1 phlebotomy armrest

AC0027 DRIP STAND
Made from chromed steel.

EXAMPLE 4
with the addition of:
• n. 2 fastening devices
• n. 1 armchair armrest
• n. 1 drip stand

AC0026 FASTENING DEVICE
Made from chromed steel (can be supplied individually)

For further information on options and accessories, please refer to the first pages.

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.



Phlebotomy Couch

This multifunctional couch is identified by code LE 311 D2W6 and it has two motors; height is adjusted electrically by means of a hand control fitted with a magnet. The backrest and leg rest sections can be adjusted via a system that is servo-assisted by gas springs and convenient levers, while the seat section is easily adjusted by means of another motor. This couch has several positions: seated, relaxation, Trendelenburg and anti-Trendelenburg. Includes wheels.

The configuration can also include the following accessories when used for phlebotomy:

AC0026 FASTENING DEVICE	n. 3 Pieces
AC0028 ARMCHAIR ARMREST	n. 1 Piece
AC0029 PHLEBOTOMY ARMREST	n. 1 Piece
AC0027 DRIP STAND	n. 1 Piece



Model shown is LE 311 D2 W6

Armchair/Couch and Transport Couch

This model has two motors; height is adjusted electrically by means of a hand control fitted with a magnet. The backrest and leg rest sections can be adjusted via a system that is servo-assisted by gas springs and convenient levers, while the seat section is easily adjusted by means of another motor. This couch has several positions: seated, relaxation, Trendelenburg and anti-Trendelenburg. Includes wheels. It is useful for patients who cannot stand up and move by themselves, who require long-term care, are bedridden for a long period of time or live in care homes, as it can help them reach the desired height to be then transferred horizontally and come to a seated position while the patient is still lying on the table, and finally be helped to move by a carer.

AC0026 FASTENING DEVICE	n. 2 Pieces
AC0028 ARMCHAIR ARMREST	n. 2 Pieces
AC0698 PUSH HANDLE	n. 1 Piece



Model shown is LE 211 D1 W4

Gynecological examination table

Our gynecological and urological examination table has three mobile sections and can have one or two motors for height adjustment and, if required, to adjust the seat section and adjust the tilt angle, with the patient lying on the table and without any strain. This examination table can be adjusted into several positions, starting from sitting at a lower height which makes it easier for the patient to climb on. The sections are adjusted by gas springs and ergonomic levers.




THESE MODELS COME WITH TWO PADDED THIGH SUPPORTS, THAT CAN BE HEIGHT-ADJUSTED AND ROTATED, AND A STAINLESS STEEL REMOVABLE BOWL AS STANDARD SUPPLY.



Model shown is LE 351 F2 W6

Test Line Multi-functional tables

Pos. 3 - 4 - 5 Type of sections



Pos. 6 Height adjustment

For models LE111, LE121

- A** Electric, with foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

For model LE311

- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control

Pos. 7 Other options


- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available



Model shown - code LE 111 A1 WT



TECHNICAL FEATURES

Number of sections	3
Max safety load (kg)	200
Height adjustment	Electric / Hydraulic
Min / max table height (cm) - LE111, LE121	60-110
Min / max table height (cm) - LE311	65-115
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+54° -52°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	82 / 82 / 89

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0037 ADDITIONAL FOOT PEDAL 1
- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

SPECIALIST ACCESSORIES:

- AC0023.W? CERVICAL CUSHION
- AC0026 FASTENING DEVICE *
- AC0027 DRIP STAND *
- AC0028 ARMCHAIR ARMREST *
- AC0029 PHLEBOTOMY ARMREST *
- AC0710.W? PADDED ARMREST *


* IMPORTANT NOTE: accessories AC0026, AC0027, AC0028, AC0029, AC0710.W? are only available for models LE111 e LE311

upholstery
color code




Pos. 3 - 4 - 5 Type of sections


211 TEST MOBIL



411 TEST MOBIL PLUS



Model shown - code LE 211 D1 WT



Pos. 6 Height adjustment

For models LE211

- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

For model LE411

- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control

Pos. 7 Other options

- I** No options

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the specific pages.



TECHNICAL FEATURES

Number of sections	3
Max safety load (kg)	200
Height adjustment	Electric / Hydraulic
Min / max table height (cm) - LE211	60-110
Min / max table height (cm) - LE411	65-115
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+54° -52°
Wheels diameter (mm)	150
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	82 / 89

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0039 ADDITIONAL HAND CONTROL 1
- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0699 PAPER ROLL HOLDER H *
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

SPECIALIST ACCESSORIES:

- AC0023.W? CERVICAL CUSHION
- AC0026 FASTENING DEVICE
- AC0027 DRIP STAND
- AC0028 ARMCHAIR ARMREST
- AC0029 PHLEBOTOMY ARMREST
- AC0710.W? PADDED ARMREST
- AC0698 PUSH HANDLE *

* IMPORTANT NOTE: "Push handle" can be combined with "Paper Roll Holder H (AC0699) only.


upholstery color code




Test Line Multi-functional tables

Pos. 3 - 4 - 5 Type of sections

141 TEST PODO



341 TEST PLUS PODO



Pos. 6 Height adjustment

For model LE141

- A** Electric, with foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

For model LE341

- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



L E | | | | | | | | | |
1 2 3 4 5 6 7 8 9

Model shown - code LE 341 F1 W6



Configuration example for a treatment table from the VISIT range.



TECHNICAL FEATURES

Number of sections	4
Max safety load (kg)	200
Height adjustment	Electric
Min / max table height (cm) - LE141	60-110
Min / max table height (cm) - LE341	65-115
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+54 -52
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	82 / 85

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

PADDED SUPPORTS AND ELEMENTS:


- AC0023.W? CERVICAL CUSHION
- AC0026 FASTENING DEVICE
- AC0027 DRIP STAND
- AC0028 ARMCHAIR ARMREST
- AC0029 PHLEBOTOMY ARMREST
- AC0710.W? PADDED ARMREST

upholstery color code




Pos. 3 - 4 - 5 Type of sections

151 TEST UROGYN



351 TEST PLUS UROGYN



Pos. 6 Height adjustment

For model LE151

- A** Electric, with foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Electric, with dual control

For model LE351

- E** Electric, two motors with foot switch
- F** Electric, two motors with hand control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation



Pos. 8 Synthetic leather design

- W** Seamless rounded edges


Pos. 9 Section upholstery color

24 options available

For further information on options and accessories, please refer to the first pages.



Accessory AC0545.W?
plug for gynecological
examination table



Model shown - code LE 351 F1 W6



Padded thigh holders and bowl are included in the standard supply.

TECHNICAL FEATURES

Number of sections	3
Max safety load (kg)	200
Height adjustment	Electric
Min / max table height (cm) - LE151	60-110
Min / max table height (cm) - LE351	65-115
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+75°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m ³)	30
Weight (kg)	82 / 92

ACCESSORIES:

ELECTRICAL ACCESSORIES:

- AC0038 ADDITIONAL FOOT PEDAL 2
- AC0040 ADDITIONAL HAND CONTROL 2
- AC0083 ADDITIONAL CONTROL 1
- AC0997 FOOT SWITCH HOLDER
- AC0578 SAFETY DEVICE
- AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

- AC0031 PAPER ROLL HOLDER F
- AC0034 STANDARD PAPER ROLL SET
- AC0035 MAXI PAPER ROLL SET
- AC0036 PROTECTIVE SHEET SET

SPECIALIST ACCESSORIES- TEST SERIES:

- AC0023.W? CERVICAL CUSHION
- AC0026 FASTENING DEVICE
- AC0027 DRIP STAND
- AC0028 ARMCHAIR ARMREST
- AC0029 PHEBOTOMY ARMREST
- AC0324 LEG SUPPORTS
- AC0710.W? PADDED ARMREST
- AC0545.W? PLUG FOR EXAMINATION TABLE


upholstery
color code



TestLine Multi-functional tables

Pos. 3 - 4 - 5 Type of sections

131 TEST ECODOP



Pos. 6 Height adjustment

- A** Electric, with foot switch
- C** Hydraulic, with double foot switch
- D** Electric, with hand control
- G** Multifunctional ergonomic control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Section upholstery color

24 options available

Configuration example for a treatment table from the VISIT range.



For further information on options and accessories, please refer to the first pages.



Model shown - code LE 131 A1 WT



TECHNICAL FEATURES

Number of sections	4
Max safety load (kg)	200
Height adjustment	Electric / Hydraulic
Minimum/maximum table height (cm)	60-110
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	54
Head section adjustment	Gas spring-assisted
Head section tilt	+75°
Wheels diameter (mm)	65
Braking system	Central
Upholstery safety rating	1 IM
Padding thickness (mm)	40
Density (kg/m³)	30
Weight (kg)	82

ACCESSORIES:

ELECTRICAL ACCESSORIES:
 AC0037 ADDITIONAL FOOT PEDAL 1
 AC0039 ADDITIONAL HAND CONTROL 1
 AC0083 ADDITIONAL CONTROL 1
 AC0997 FOOT SWITCH HOLDER
 AC0017 AUXILIARY BATTERY

FOR PROTECTION AND HYGIENE:

AC0031 PAPER ROLL HOLDER F
 AC0034 STANDARD PAPER ROLL SET
 AC0035 MAXI PAPER ROLL SET
 AC0036 PROTECTIVE SHEET SET

SPECIALIST ACCESSORIES- TEST SERIES:

AC0023.W? CERVICAL CUSHION
 AC0545.W? PLUG FOR EXAMINATION TABLE

upholstery color code



AluX Eva



With this range, Chinesport is branching out into beauty, spa and massage tables for the first time, and introduces some additional optional features such as the sections heated at three different temperatures and the synthetic leather upholstery that covers the sides to produce a softer padding effect on the contact surface. Among other noteworthy options, there is a possibility to choose a led lights configuration for combining chromotherapy effects. This is why this model comes with more identification codes than any other height-adjustable table, and even when there is only one option of section width or padding density/thickness clients can request personalized features and options. These tables have two motors, so they can achieve the Trendelenburg and anti-Trendelenburg position via a dual telescopic column system.

POSITION 4 - Adjusting the head section



1 **Head section without armrests, adjustable by gas spring** with lever control. This section can be tilted up to + 70°.



2 **Head section without armrests, adjustable by electric motor** with hand control with magnet. This section can be tilted up to + 70°.



3 **Head section with armrests, adjustable by gas spring** with lever control. This section can be tilted up to + 70°.



4 **Head section with armrests, adjustable by electric motor** with hand control with magnet. This section can be tilted up to + 70°.


POSITION 5 - Adjusting the legs-trunk sections



1 **Combined trunk-leg section.** Only one section for trunk and legs which can be adjusted by an electric motor in positive or negative tilt +/- 12°.



3 **Movable sections, adjustable by gas spring** - Two movable sections obtaining a relax position, they are adjustable by gas strut with lever control. It is possible to adjust the tilt of all sections together + / - 12°.



4 **Movable sections, adjustable by electric motor** - Two movable sections obtaining a relax position, they are adjustable by electric motor with magnet. It is possible to adjust the tilt of all sections together + / - 12°.

POSITION 6 - Adjusting the frame in height



F **With hand control**
The adjustment is done by means of a hand control with magnet;

M **With hand control on both sides**
This option includes not only the magnet control by hand, but also a fix control on both sides.

POSITION 7 - Transferring wheels and led illumination



1 Without wheels, without led illumination



2 With transferring wheels

This table can be provided with retractable wheels with centralized brake system lever activation. The wheels enable easy transferring without person lying on it.



8 With led illumination

9 With wheels, with led illumination

The ALUX therapy table is available with led illumination with a possible combination of 640.000 colours. Intensity can be programmed up to 20 ways of illumination by means of a specific wireless control. These lights can provide chromotherapy effects, and they help make the table become harmonically integrated in the place where it is.

POSITION 8, 9, 10, 11 - Width, density and thickness



075 Standard width cm 75

3 Thickness mm 50, Density 30 kg/m³

The width of the sections is of cm 75 standard, while the padding is mm 50 minimum with a 30 kg / m³ density. Our upholstering center can make personalizations, after verifying the feasibility. In such cases an adjustment of the price would be done.

POSITION 12 - Heating



A Padded sections without heating

B Padded sections with heating

For a greater comfort during treatment it is possible to set the bed with a heating function of the sections. Thus, the sections will always be warm and pleasant at touch for a more relaxing experience.

POSITION 13 - Breathing hole



1 Head section without breathing hole



2 Head section with breathing hole

POSITION 14, 15 - Synthetic leather finishing and upholstery colours



X Sections with side band

The sections are made with a side band which makes padding softer. The type of synthetic leather is Giotto type collection. 24 colours of upholstery are available.

Giotto range

colour code upholstery



Pos. 4 Head section

- 1** without armrests, gas spring adjustment
- 2** without armrests, electric adjustment
- 3** with armrests, gas spring adjustment
- 4** with armrests, electric adjustment

Pos. 5 Trunk and legs sections

- 1** One section, tilt only
- 3** Sections adjustable by gas springs
- 4** Sections adjustable electrically

Pos. 6 Height adjustment

- F** Electric, two motors with hand control
- M** Electric with hand control and keypad

Pos. 7 Wheels and led lights

- 1** without wheels, without led lights
- 2** with wheels
- 8** with led lights
- 9** with wheels and led lights

Pos. 8 - 9 - 10 Section width

- 0 7 5** width 75 cm

Pos. 11 Density / Thickness

- 3** Density 30 kg/m³ Height 50 mm

Pos. 12 Heating

- A** Not heating
- B** With heating

Pos. 13 Nose / mouth breathing hole

- 1** Without breathing hole
- 2** With breathing hole

Pos. 14 Type of synthetic leather design


- X** Squared edges

Pos. 15 Section upholstery color

24 options available

For further information on options and accessories, please refer to the specific pages.

L	A	2								0	7	5	3						
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15					



Model shown - code LA2 11M1 0753 A1 XR



TECHNICAL FEATURES

Number of sections	2/4/6
Max safety load (kg)	180
Height adjustment (dual telescopic column)	Electric
Minimum/maximum table height (cm)	60-100
Mains supply	230V 50/60 Hz / 24V
Head section adjustment	Gas spring / Electric
Head section tilt	+70°
Trunk-legs sections adjustment	Gas spring / Electric
Trendelenburg / Counter-trendelenburg	Electric
Wheels with central braking system	Optional
Upholstery safety rating	1 IM
Padding thickness (mm) / Density (kg/m) - X	50 / 30
Weight (kg)	85



Integrated plugs for power supply.

ACCESSORIES:
AC0034 STANDARD PAPER ROLL SET
AC0035 MAXI PAPER ROLL SET
AC0784 PROTECTIVE SHEET SET
AC0287 CLASSIC PAPER ROLL HOLDER
AC0770.3X? BREATHING HOLE PLUG
AC0024.W? FACE CUSHION

upholstery color code




Configuration example n. 1

L A 2 1 1 M 1 0 7 5 3 A 1 X R
 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15



1 Pos. 13 - Head section without breathing hole



A Pos. 12 - Padding without heating



1 Pos. 4 - Head section without armrests, gas spring adjustment



1 Pos. 5 - One trunk and legs section, +/-12° tilt adjustment





M Pos. 6 - Electric height adjustment with double control



1 Pos. 7 - Without wheels, and without led lights

Configuration example n. 2

L A 2 4 4 M 9 0 7 5 3 B 2 X A
 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15




2 Pos. 13 - Head section with breathing hole



B Pos. 12 - Padding with heating



4 Pos. 4 - Head section with armrests, electric adjustment



4 Pos. 5 - Trunk and legs sections electrically adjustable



M Pos. 6 - Electric height adjustment with double control



9 Pos. 7 - With wheels and led lights

upholstery color code



FOR PROTECTION AND COMFORT



AC0034 STANDARD PAPER ROLL SET

The set includes 6 x 60 cm wide rolls. Each roll is 85 meters long.



AC0035 MAXI PAPER ROLL SET

The set includes 2 x 90 cm wide rolls. Each roll is 100 meters long.



AC0784 PROTECTIVE SHEET SET

This set contains 10 synthetic fabric sheets with elasticated border. This sheet can be used only for Alux and Sinthesi series.



AC0287 CLASSIC PAPER ROLL HOLDER

This holder can only be attached to the head section of ALUX and MINITOP models and can only be used with paper rolls that are up to 61 cm wide.



AC0770.3X ? BREATHING HOLE PLUG

Dimensions: L 31 x P 23 x H 6



AC0024.W ? FACE CUSHION

Dimensions: L 31 x P 23 x H 6



XCAT0020131 UPHOLSTERY COLOR BOOK

Upon request from design companies or individual professionals who would like to check the actual synthetic leather color before making a purchase, we are happy to send our samples.



An overview on the main features of this new Chinesport special table for aesthetic treatments.


? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available



Basic tables



Examination or treatment tables with a sturdy metal or wooden frame; dimensions and number of segments vary but typically have fixed height. However, some models can be height-adjusted manually to pre-determined positions, while others can be folded. The section upholstery has rounded edges, which allows patients to come closer, while the head section can be ordered in the version with nose & mouth breathing hole. These tables can be used in doctor's surgeries, for example. The MINITOP range is a recent addition to this group and is a basic model for examinations only, with two sections that can be adjusted electrically and a very simplified design.



Accessories for basic tables

FOR PROTECTION AND HYGIENE



AC0031 PAPER ROLL HOLDER F

A metal support with adjustable width from 50 cm to 65 cm that is attached to the upper part of a couch. Recommended for use with paper rolls from set AC0034. This holder can be attached to basic tables on either side, head or foot.



AC0287 CLASSIC PAPER ROLL HOLDER

This holder can only be attached to the head section of ALUX and MINITOP models and can only be used with paper rolls that are up to 61 cm wide.



Y740 ROLL HOLDER

This metal frame holds the paper rolls used on examination, massage and treatment tables and can be easily moved.



AC0034 STANDARD PAPER ROLL SET

The set includes 6 x 60 cm wide rolls. Each roll is 85 meters long.



AC0035 MAXI PAPER ROLL SET

The set includes 2 x 90 cm wide rolls. Each roll is 100 meters long.



AC0036 PROTECTIVE SHEET SET

This set contains 10 synthetic fabric sheets with elasticated border. This sheet can be used on any examination and treatment table. Not compatible with couches larger than cm 100 or more.

SUPPORTS AND UPHOLSTERY



AC0020.W NOSE/MOUTH BREATHING HOLE PLUG

For treatment tables that provide for the option of a standard nose/mouth breathing hole.



12269 HEADREST

Tilt-adjustable cervical cushion that can be fitted on all basic treatment tables using a special frame. Available only in gray. Not compatible with couches larger than cm 100 or more.




12270 ARMREST

This support helps to release the load from arms and can be ordered together with the headrest code 12269. Available only in gray. Not compatible with couches larger than cm 100 or more.



Pos. 3 - 4 - 5 Type of table top

E10 MINITOP B



Pos. 6 Height adjustment

- A** Electric, with foot switch
- D** Electric, with hand control

Pos. 7 Other options

- 1** No options
- 2** With wheels for transportation

Pos. 8 Synthetic leather design

- W** Seamless rounded edges

Pos. 9 Upholstery color

24 options available




Cushion cod. AC0025, improves the comfort when in prone position (see page 10 for more info)



Detail of the ergonomic control for the head section adjustment.

For further information on options and accessories, please refer to the first pages.



Model shown LVE10 D1 W3



TECHNICAL FEATURES

Number of sections	2
Max safety load (kg)	150
Height adjustment	Electric
Minimum/maximum table height (cm)	60-90
Mains supply	220V 50-60Hz / 24V 70VA
Actuator thrust (N)	6000
IP Protection	66
Head section adjustment	Gas spring-assisted
Head section tilt	0° +75°
Wheels diameter (mm)	5
Braking system	Not applicable
Upholstery safety rating	1 IM
Padding thickness (mm)	30
Density (kg/m³)	30
Weight (kg)	50



Detail of the wheel for easy transfer

ACCESSORIES

- AC0034 2-PLY PAPER ROLL
- AC0031 PAPER ROLL HOLDER
- AC0287 CLASSIC PAPER ROLL HOLDER

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available



Basic tables *metal frame*




01250.W?

METAL 1

An examination and treatment table with two sections and fixed height at 78 cm. The standard section width is 68 cm, while the head section can have manual positive tilt adjustment up to +45°. The head section in this model has a nose/mouth breathing hole.


A breathing hole plug can be ordered separately as an accessory. Max safety load: 200 kg; Weight: 40 kg



01650.W?

METAL BASIC 1

This model is similar to item code 01250.W, the only difference being that the head section does not have a nose/mouth breathing hole.




01252.W?

METAL L

An examination and treatment table with two sections and fixed height at 78 cm. Sections have a larger width of 80 cm, while the head section can have manual positive tilt adjustment up to +45°. The head section in this model has a nose/mouth breathing hole. A breathing hole plug can be ordered separately as an accessory.


Max safety load: 200 kg; Weight: 43 kg



01652.W?

METAL BASIC L

This model is similar to item code 01252.W, the only difference being that the head section does not have a nose/mouth breathing hole.




01251.W?

METAL 3

An examination and treatment table with four sections and fixed height at 78 cm. The standard section width is 68 cm, while the mobile sections can be adjusted manually. The head section in this model has a nose/mouth breathing hole.


A breathing hole plug can be ordered separately as an accessory. Max safety load: 200 kg; Weight: 43 kg



01651.W?

METAL BASIC 3

This model is similar to item code 01251.W, the only difference being that the head section does not have a nose/mouth breathing hole.



Please refer to the introduction page of this section for further details on our accessories.

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.






01256.W ?

METAL H100

An examination and treatment table with one single section and fixed height at 78 cm. Larger dimensions for a total surface area that is 100 x 200 cm.

Max safety load: 250 kg; Weight: 58 kg




01257.W ?

METAL HT100

An examination and treatment table with two sections and fixed height at 78 cm. Larger dimensions for a total surface area that is 100 x 200 cm. The head section can have manual positive tilt adjustment up to +55°.

Max safety load: 250 kg;

Weight: 60 kg




01253.W ?

METAL L100

An examination and treatment table with one single section and fixed height at 50 cm. Larger dimensions for a total surface area that is 100 x 200 cm.

Max safety load: 250 kg;

Weight: 53 kg




01254.W ?

METAL LT100

An examination and treatment table with two sections and fixed height at 50 cm. Larger dimensions for a total surface area that is 100 x 200 cm.

The head section can have manual positive tilt adjustment up to +55°. Max safety load: 250 kg; Weight: 55 kg



01423.W ?


METAL L200

Metal examination and treatment table with two connected sections, for a total surface area that is 200 x 200 cm.

Fixed table height at 50 cm.

Max safety load: 250 kg;

Weight: 106 kg



Please refer to the introduction page of this section for further details on our accessories.

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.




Basic tables *wooden frame*



01258.W?

RAMIN 1


Wooden examination and treatment table with two sections and fixed height at 78 cm. The standard section width is 68 cm, while the head section can have manual positive tilt adjustment up to +45°. The head section in this model has a nose/mouth breathing hole. A breathing hole plug can be ordered separately as an accessory. Max safety load: 200 kg; Weight: 40 kg



01658.W?

RAMIN BASIC 1


This model is similar to item code 01258.W, the only difference being that the head section does not have a nose/mouth breathing hole.



01260.W?

RAMIN L


Wooden examination and treatment table with two sections and fixed height at 78 cm. The standard section width is 80 cm, while the head section can have manual positive tilt adjustment up to +45°. The head section in this model has a nose/mouth breathing hole. A breathing hole plug can be ordered separately as an accessory. Max safety load: 200 kg; Weight: 45 kg



01660.W?

RAMIN BASIC L


This model is similar to item code 01260.W, the only difference being that the head section does not have a nose/mouth breathing hole.



01259.W?

RAMIN 3


Wooden examination and treatment table with four sections and fixed height at 78 cm. The standard section width is 68 cm, while the mobile sections can be adjusted manually. The head section in this model has a nose/mouth breathing hole. A breathing hole plug can be ordered separately as an accessory. Max safety load: 200 kg; Weight: 43 kg



01659.W?

RAMIN BASIC 3

This model is similar to item code 01259.W, the only difference being that the head section does not have a nose/mouth breathing hole.



Please refer to the introduction page of this section for further details on our accessories.

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.






01263.W ?

RAMIN H100


Wooden examination and treatment table with one single section and fixed height at 78 cm. Larger dimensions for a total surface area that is 100 x 200 cm. Max safety load: 250 kg; Weight: 58 kg



01264.W ?

RAMIN HT100


Wooden examination and treatment table with two sections and fixed height at 78 cm. Larger dimensions for a total surface area that is 100 x 200 cm. The head section can have manual positive tilt adjustment up to +55°. Max safety load: 250 kg; Weight: 60 kg



01261.W ?

RAMIN L100


Wooden examination and treatment table with one single section and fixed height 50 cm. Larger dimensions for a total surface area that is 100 x 200 cm. Max safety load: 250 kg; Weight: 53 kg



01262.W ?

RAMIN LT100

Wooden examination and treatment table with two sections and fixed height at 50 cm. Larger dimensions for a total surface area that is 100 x 200 cm. The head section can have manual positive tilt adjustment up to +55°. Max safety load: 250 kg; Weight: 55 kg



Please refer to the introduction page of this section for further details on our accessories.

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.




Basic tables *wooden frame*



01265.W ?

VOJTA


Wooden examination and treatment table with one mobile section, head section with manual adjustment system. The head section has a nose/mouth breathing hole. With a convenient storage shelving underneath. Fixed height 78 cm
Max safety load: 200 kg; Weight: 50 kg



01653.W ?

VOJTA BASIC


This model is similar to item code 01265.W, the only difference being that the head section does not have a nose/mouth breathing hole.



02988.W ?

VOJTA 1R


Height-adjustable, wooden examination and treatment table with one mobile section. Head section with manual adjustment system and nose/mouth breathing hole with plug. With a convenient storage shelving underneath. Min/max height 63 ÷ 87 cm.
Max safety load: 200 kg; Weight: 47 kg



01654.W ?

VOJTA BASIC 1R


This model is similar to item code 02988.W, the only difference being that the head section does not have a nose/mouth breathing hole.



02989.W ?

VOJTA 3R


Height-adjustable, wooden examination and treatment table with three mobile sections. Head and leg sections with manual adjustment system. The head section has a nose/mouth breathing hole with plug. With a convenient storage shelving underneath. Min/max height 63 ÷ 87 cm.
Safety load 200 kg; Weight: 52 kg



01655.W ?

VOJTA BASIC 3R

This model is similar to item code 02989.W, the only difference being that the head section does not have a nose/mouth breathing hole.



Please refer to the introduction page of this section for further details on our accessories.

? CHOICE OF UPHOLSTERY COLOR - important: always specify the upholstery code along with the chosen item when this option is available.






11975.W2

CHIROPRACTIC TABLE


Chiropractic treatment table with folding, height-adjustable leg supports. Min/max height 61 - 68 - 74 cm. The elbow support section can be extended. Washable, synthetic leather upholstery, only available in blue W2. Max safety load: 150 kg. Weight: 34 kg



11970.W2

C.E.R.D.O. THERAPY TABLE


Treatment table with folding, height-adjustable leg supports. Min/max height 60 ÷ 80 cm. High density foam padded top upholstered with washable synthetic leather, only available in blue W2. Max safety load: 150 kg. Weight: 34 kg



11861.?

PORTABLE THERAPY TABLE

Aluminum treatment table, with two sections, folds into a suitcase. Fixed height 75 cm. This treatment table is available with black or blue synthetic leather upholstery. Safety load 170 kg; Weight: 16 kg



11863.?

PORTABLE THERAPY TABLE LARGE


Same features of model code 11861, but with width 70 cm. Weight 19 kg.



11860.?

PORTABLE THERAPY TABLE A

Aluminum treatment table, with three sections, folds into a suitcase. Fixed height 75 cm. Headrest with manual adjustment system (+ 80°) and nose/mouth breathing hole. This treatment table is available with black or blue synthetic leather upholstery. Safety load 170 kg; Weight: 17 kg



11862.?

PORTABLE THERAPY TABLE A LARGE

Same features of model code 11860, but with width 70 cm. Weight 19 kg.



ACCESSORIES

11864 TRANSPORT BAG 56

Transport protection cover suitable for foldable therapy table featured by width 56 cm.

11867 TRANSPORT BAG 70

Transport protection cover suitable for foldable therapy table featured by width 70 cm.



Basic tables

CUSHIONS FOR POSTURE


Postural cushions are available in various shapes and sizes to meet the professional's varied requirements to position a patient correctly. Our cushions are made from variable density foam and have a synthetic leather cover which is fire-resistant and anti-bacterial, easy to clean and durable. Standard density is 21 kg/m³. Cushion covers have a zip for easy removal, and come in a choice of upholstery from an assortment of 24 color options.



Please refer to the first pages of this section for further details.

Other furniture

We have a range of padded stools in various shapes that can be chosen by the professional among several options. Our stools can also be made in the same upholstery color as the Chinesport therapy tables, so there are 24 colors to choose from. Below is an interior design solution for medical and physical therapy practices including a series of items and accessories such as wooden stools, chairs, office chairs, screens, medical equipment trolleys. The step plates with one or two steps make it easier for patients to be transferred onto fixed height therapy tables.



CUSTOM STOOLS

Series of swivel stools with five swivel wheels and 60 cm diameter base. Height is adjusted by means of a gas spring mechanism and lever. Models come with two different height adjustment mechanisms. Our stools are manufactured using materials that ensure excellent seat stability over time. Bushings are also used to attach the seat to the frame, which allows for parts to be easily replaced when worn out or if the padding gets damaged, or even if you just want a new upholstery color. A wide choice of upholstery colors is available. 3-year guarantee.

Pos. 4 Synthetic leather design



1 Seamless rounded edges



2 Squared edges

Pos. 5 Backrest



1 Without backrest



2 With backrest

Pos. 6 Height adjustment



L from 46 to 60 cm



H from 60 to 90 cm

Pos. 7 Type of frame



1 Black plastic



2 Chromed



MADE-TO-MEASURE SEAT

We can personalize the seat height, density and shape. Feel free to send us your request by e-mail to our address: export@chinesport.it

Pos. 3 Basic model

- 1 Miguel
- 2 Ines
- 3 Saddle stool
- 4 Ronaldo

Pos. 4 Synthetic leather design

- 1 Seamless rounded edges
- 2 Squared edges

Pos. 5 Backrest

- 1 Without backrest
- 2 With backrest

Pos. 6 Height adjustment

- L from 46 ÷ 60 cm
- H from 60 ÷ 90 cm

Pos. 7 Type of frame

- 1 Black plastic
- 2 Chromed

Pos. 8 Type of synthetic leather

- W Giotto Collection

Pos. 9 Upholstery color

24 options available

For further information on upholstery characteristics, please refer to the first pages.



Model shown - code CS1 12L2 WE

MIGUEL

C	S	1						W
1	2	3	4	5	6	7	8	9



Model shown - code CS2 21H1 W3

INES

C	S	2	2					W
1	2	3	4	5	6	7	8	9



Model shown - code CS3 11H1 WT

SADDLE STOOL

C	S	3	1	1				W
1	2	3	4	5	6	7	8	9



Model shown - code CS4 11L1 WQ

RONALDO

C	S	4	1	1				W
1	2	3	4	5	6	7	8	9

upholstery color code



Other furniture



11815

FIXED STOOL

Chromed steel with beechwood seat.
SEAT Ø 35 X 48h cm



11786

SWIVEL STOOL 2

Chromed steel with beechwood seat and 5 swivel wheels. Height adjustable with screw mechanism.
60 x 40 ÷ 53h cm



XEC001

RONDO' STOOL

A stool on a chromed base fitted with 5 casters. The backrest is designed for use as armrest as well. Height adjustable 62 - 82 cm.



02077

CHAIR WITHOUT ARMRESTS

Natural painted beechwood.
52 x 50 x 83h cm
Weight 7 kg



02076

CHAIR WITH ARMRESTS

Natural painted beechwood.
52 x 50 x 83h cm
Weight 9 kg



12880

WOODEN STOOL

50 x 36 x 47h cm



12855

WOODEN FOOTREST

Solid beechwood.
32 X 22 X 16h cm



11810

SQUARE STOOL

Made from section tube, which offers maximum stability.
40 x 40 x 45h cm



XP851

STEP PLATE 1

Made from chromed steel tube.
1 plastic step, 38 x 20 cm.
25 x 38 x 23h cm
Capacity 150 kg



XP852

STEP PLATE 2

Made from chromed steel tube.
2 plastic steps, 38 x 20 cm.
50 x 38 x 44h cm
Capacity 150 kg



34630

WAITING ROOM CHAIR 1

Frame made from coated, oval steel; padded seat and backrest but no armrests.
58 x 53 x 80h cm



01283

WAITING ROOM CHAIR 2

Frame made from coated, oval steel; padded seat and backrest with removable armrests.
58 x 64 x 80h cm



34650

3-SEAT CHAIR

Side-by-side frame made from coated, oval steel; padded seat and backrest.
160 x 60 x 80h cm



34620

OFFICE CHAIR WITHOUT ARMRESTS

Height-adjustable with swivel base and 5 swivel wheels. Adjustable backrest that can be used to rest on at all times.
Base Ø 60 x 44÷57h cm



34610

OFFICE CHAIR WITH ARMRESTS

Height-adjustable with swivel base and 5 swivel wheels. Adjustable backrest that can be used to rest on at all times.
Base Ø 60 x 44÷57h cm



Y766.P53

VERONICA

Ergonomic chair complete with knee support and adjustable armrests. Tilting backrest. Gas lift mechanism.
Base Ø 60 x 53÷66h cm

Upholstery color

P53

Other furniture



02112
OLMO DESK
80 x 140 x 73h cm



02113
OLMO CHEST OF
DRAWERS
Three drawers with lock, wheels
and metal handles.
45 x 48 x 67h cm



XM2520
TROLLEY 2
Aluminium tube frame that can be
disassembled. Fitted with 4 swivel
wheels. 2 white bilaminat shelvings
with rounded edges. 65 x 37 x 83h
cm



XM2530
TROLLEY 3
Aluminium tube frame that can be
disassembled. Fitted with 4 swivel
wheels. 3 white bilaminat shelvings
with rounded edges. 65 x 37 x 90h
cm




34030
RUBBISH BIN
12-liter container,
plastic on the inside.
Ø 25 x 38h cm



02045
SPIRAL BOX 2
Wall-mounted paper roll
dispenser. Only for tear-off
paper.



13190
PAPER 1
2-ply pure cellulose, tear-off,
180 m. Pack contains 6 rolls.
For "Spiral Box 2" code 02045.



XM16506

MOBILE SCREEN

Screen with 6 folding laminate panels plus one support panel. The wheels on the base allow for it to be moved to the desired location as required. Aluminum frame. Braking casters. Screen useful length 180 cm. Overall height 165 cm.

XM18506

WALL-MOUNTED SCREEN

Screen with 6 folding laminate panels that attaches to the wall using a special mounting rail. Aluminum frame. Braking wheels. Screen useful length 180 cm. Overall height 185 cm.

XM16508

MOBILE SCREEN

Same as item XM16506 but with 8 folding panels, reaching a useful length of 240 cm.

XM18509

WALL-MOUNTED SCREEN

Same as item XM18506 but with 8 folding panels, reaching a useful length of 240 cm.



34070

METAL COAT HANGER

Made from coated steel tube. With umbrella stand. Base 58 x 175h cm

02103

3-SECTION SCREEN

Chromed steel tube frame. Ø 2.2 cm with PVC. 3 sections x 50 cm x 164h



Other furniture *Medical equipment trolleys*









XM2703 3-DRAWER TROLLEY

Storage trolley with an aluminum and ABS frame. It has 2 drawers that are 8 cm tall, 1 drawer that is 16 cm tall, and one open shelving. Four swivel wheels, two of which with braking mechanism. Central locking system with key. Dimensions 76x50x87 h cm. NOTE: accessories represented in picture are available for purchase separately.




XM2704 4-DRAWER TROLLEY


Storage trolley with an aluminum and ABS frame. It has 2 drawers that are 8 cm tall, 1 drawer that is 16 cm tall, and 1 drawer that is 24 cm tall. Four swivel wheels, two of which with braking mechanism. Central locking system with key. Dimensions 76x50x87 h cm. NOTE: accessories represented in picture are available for purchase separately.




Examples of accessories. To be attached on the short sides of a trolley using a height-adjustable bar



XM0004
SHARPS CONTAINER



XM0005
CLOSED RUBBISH BIN



XM0006
GLOVES DISPENSER

ACCESSORIES

- XM0004** SHARPS CONTAINER
- XM0005** CLOSED RUBBISH BIN
- XM0006** GLOVES DISPENSER
- XM0007** BOTTLES HOLDER
- XM0008** MULTIFUNCTIONAL STORAGE
- XM0009** DOCUMENT HOLDER
- XM0010** MEDICAL RECORDS HOLDER
- XM0011** OPEN RUBBISH BIN
- XM0012** CATHETER HOLDER
- XM0013** SUCTION TUBE HOLDER