




Your way to
fitting excellence

The audiometer

655g of concentrated knowledge



Unique design: Compact and designed to support your needs.



Connectivity: The vast connectivity possibilities make it the natural central point of the fitting solution.

Under desk mounting: Mounting the audiometer with the built-in Hi-Pro 2 under the desk helps tidy up cabling and provides convenient access to the Hi-Pro plugs.


Aurical® Aud is the world's only true fitting audiometer. It integrates an audiometer, counseling tools, power amplifier, USB hub and optional Hi-Pro 2 with High Speed Mode in a modern, ultra-compact design. It's a complete, portable audiometer that offers outstanding flexibility. It is designed to fit into all kinds of environments, inside and outside your clinic.

The audiometer is easy to use too. The exceptionally intuitive Otosuite® interface has an on-screen control panel and easy-to-access User Test functions that can be customized to suit individual needs. It ensures fast operation and consistently accurate results.


High quality sound is important for your fitting results. The unique RoomTune™ functionality turns Aurical Aud into a compensated sound system allowing you to present more accurate and reliable sound from your fitting modules.

THE FACTS

- AC, BC and SF audiometry (built-in power amplifier)
- Tone audiometry testing up to 12.5 kHz, including FRESH noise™
- Loudness Scaling (optional)
- Integrated speech material
- SAL, TEN, QuickSIN tests (some optional)
- Masking assistant
- Sound adjustment with RoomTune™
- Built-in USB Hub
- Optional Hi-Pro 2 w/ High Speed Mode
- Otosuite with User Tests
- Intuitive on-screen control panel
- One-click data logging
- Selection of report templates



Various view options allow you to choose what is best for you and your patient. Go from simple to advanced view in 2 mouse clicks.



Quick Select

AC: R B L

Mask:

BC:

Mask:

SF:

Mask:

Condition: Threshold MCL UCL

Test Options

dB Step: 1 5

AC Transducer: Insert Phone

Control: Stim. On Stim. Lock Tracking Pulsed

Stimulus: Tone Warble Stimulus 2

Monitor / Level

MONITOR: Stim Mask Talk back

LEVEL: Mic Source A Source B

Quick Select

Test Ear: R B L

Measurement

AC BC SF

SDT SRT MCL UCL WRS/SRS

Unaided Masking On

Test Options

dB Step: 1 5

AC Transducer: Insert Phone

Control: Stim. On Stim. Lock Tracking

Input: Int. File Source A

Monitor / Level

MONITOR: Stim Mask Talk back

LEVEL: Mic Source A Source B

The on-screen control panels in Aurical Aud and Otosuite are intuitive and fast.

Test Selector

User Test:

- Adult - Extensive Aud
- Adult - Initial Aud
- HI Benefit Aided Speech SF
- HI Benefit Aided Speech SF in Noise
- Pediatric - Extensive Aud
- Pediatric - Initial Aud

Special Test:

- 1. Full Data Overview
- 2. ABLB
- 3. SISI
- 4. Weber (+ Rinne Data Entry)
- 5. Stenger
- 6. Tone Decay
- Default USA

Cancel Help

Otosuite® lets you administer your preferences easily. Make the settings once and save as User Tests for later. Share your User Tests with your colleagues and make Otosuite ensure a smooth and efficient way of working.

The Aurical family

Fully integrated with Otosuite



The Otosuite® software integrates all the elements of Aurical®. It's intuitive, easy to learn and easy to use. It supports the entire client journey, from diagnostics and counseling, to fitting and verification. Simple to use user test functionality ensures configurable, consistent and smooth workflows that can be customized to suit individual needs. And you work in a single software environment across modules. The result is a better fitting experience – for you and for your client.

See the videos about the complete fitting solution at hearing-balance.natus.com/aurical

About Natus

A leader in hearing and balance products and solutions, Otometrics was acquired in 2017, and is now known as Natus. For almost 60 years, we have been helping clinicians improve the quality of life for their clients and patients by delivering expert knowledge, reliable solutions and services and trusted partnerships. We continue to develop, manufacture and market audiological, otoneurologic and vestibular instrumentation in more than 80 countries under trusted brand names, including Madsen®, Aurical®, ICS®, Otoscan® and Bio-logic®. Natus Medical Incorporated is the global brand clinicians depend on when diagnosing and treating disorders of the central nervous and sensory systems to improve patient outcomes. For more information, please visit natus.com.

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 7-26-8401-EN/02. P/N 7-26-84001-EN.

natus®

Natus Medical Incorporated

natus.com