



## Aurical® Otocam 300



### **Effective counseling and diagnostics**


For truly successful counseling, the key is to involve the patient through each step of the journey by helping them understand the procedure, diagnosis, treatment options and fitting process. Aurical Otocam 300 gives you the perfect means to start building the bridge between you and your patient.

Otocam is a tool that was designed for ease-of use, allowing you to better counsel your patient. The powerful light source lets you present detailed images directly to your computer or display. The advanced, yet intuitive, software allows the user to do a quick visual inspection for trouble-shooting or documentation throughout the process. With the built-in heating system, you will experience reduced delays attributed by condensation. It's the perfect diagnostic tool to help maximize your workflow efficiency.

There are many reasons why you should be using the Aurical Otocam 300. With easy installation and use, it's one tool that will change the way you provide diagnostic care to your patients.

### FOR DIAGNOSTICS AND COUNSELING

- Unique ease-of-use
- Unique built-in heater to minimize condensation
- No rechargeable batteries
- Powerful light source for crisp images
- Lightweight and portable design
- Intuitive software
- Plug and play installation
- Specula with and without cerumen management
- Noah compatible through Otosuite®
- No special drivers required for broad compatibility



Aurical Otocam 300 can be accessed directly from Aurical Freefit for probe tube placement and quick visual inspections during PMM.



Advanced software provides detailed, full-color videos and imagery on your computer or display.



A unique integration of the video otoscope and software allows you to remote-control the most important functions.

# The Aurical family



The Otosuite® software integrates all the elements of Aurical®. It's intuitive, easy-to-learn and easy-to-use. It supports the entire patient journey, from diagnostics and counseling to fitting and verification. User-friendly functionality ensures configurable, consistent and smooth workflows that can be customized to suit individual needs, and you work in a single software environment across modules. The result is a better fitting experience — for you and for your patient.

## See the videos about the complete fitting solution at hearing-balance.natus.com/aurical

#### **About Natus**

A leader in hearing and balance products and solutions, Otometrics was acquired in 2017, and is now known as Natus. For almost 60 years, we have been helping clinicians improve the quality of life for their clients and patients by delivering expert knowledge, reliable solutions and services and trusted partnerships. We continue to develop, manufacture and market audiological, otoneurologic and vestibular instrumentation in more than 80 countries under trusted brand names, including Madsen®, Aurical®, ICS®, Otoscan® and Bio-logic®. Natus Medical Incorporated is the global brand clinicians depend on when diagnosing and treating disorders of the central nervous and sensory systems to improve patient outcomes. For more information, please visit natus.com.

### Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. P/N 7-26-90101-EN/02.

