Robust and reliable treadmill for sports medicine











Highlights

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

Unique design

The treadmill is designed as part of the Lode design line to fit perfectly in the environment it is used.

Various handrails available

Various handrails are optional available, making the treadmill suitable for your specific stress test setting.

Smooth acceleration

All treadmills have a smooth acceleration from 0 km/h to the start speed.

Extreme performance option available

Various options make this treadmill the gold standard for top athletes:

- high speeds
- negative elevation
- fallstop





Robust and reliable treadmill for sports medicine



Modern designed treadmill specifically designed for sports medical applications. This treadmill can be connected to external ECG- and pulmonary devices as well as to a Lode software network (LEM or LCRM). The Valiant 2 sport offers a smooth acceleration from 0 km/h and is continuously adjustable in a range of 1 – 25 km/h (0.6 – 15.5 mph). Thanks to the unique low design of the cover plate it is possible to place a mirror or camera in front of the treadmill which makes it easier to monitor the gait of the test subject. The Valiant 2 sport has a standard running surface of 70 x 190 cm, 0 - 25% elevation and a full front handrail. The treadmill is equipped with a 7" control unit with Touchscreen for manual control. Other dimensions are available upon request.

A USB A-B cable only for service purposes will be standard delivered with the product.

To connect LEM or LCRM you need a special interface cable that can be ordered under part number 930930.

Features



Low noise

Material choices, refined components and accurate manufacturing techniques lead to low noise.



Versatile Interfacing

Various interface protocols guarantee perfect communication with all commonly known stress ECG and spirometry equipment

17 cm

Extreme low step-up height

To allow people to safely and comfortably step up the treadmill it is important to have a very low step-up height

0.1 km/h

Small adjustment steps

The speed of the treadmill can be adjusted in the smallest steps of only 0.1 km/h!



Low cover plate

The motor compartment of the treadmill is designed in such a way that the cover is only marginally higher than the belt surface. This allows for low camera positions in case the treadmill is used for gait analysis purposes. Also a therapist has the best possible view on feet and lower extremities.



Robust design

The product is designed to withstand continuous heavy use by subjects in most weights and sizes.



Downhill walking as an option (-10%)

This treadmill can be executed with 10% negative elevation. This allows for downhill walking which is extremely useful for rehabilitation of certain injuries.



Optional integrated BPM and SpO2

The treadmill can be extended with a stable and reliable blood pressure module and SpO2 measurement.



Robust and reliable treadmill for sports medicine

Valiant 2 sport can a.o be extended with the following options:



USB to Serial converter

Easy connection



Partnumber: 226012

Universal treadmill Arm Support

Comfort for both test subject and test



Partnumber: 945805

Emergency Stop Button

Ultimate safety



Partnumber: 945804

Bottle holder

Convenient beverage storage



Partnumber: 401200

Extension for emergency lanyard

Original Accessory



Partnumber: 945931

Autospeed

Flexibility in exercising



Partnumber: 945840

SpO2 for control unit with touch panel (extra long cable)

Oxygen saturation



Partnumber: 945822

Reverse Walking for Valiant 2

Simple switching between forward and



Partnumber: 938842

Handrail, front singlesided Valiant 2

Obstacle free observation of the test



Partnumber: 938826

Handrail, front singlesided - Pediatric -Valiant 2

Obstacle free observation of the test



Partnumber: 938816

Safety Belt and Fall stop for Valiant 2 - 70 cm

Extra safety



Partnumber: 938823

Negative elevation -10% for Valiant 2 - 70 cm

Downhill running



Partnumber: 938825

Handrail, Front -Pediatric for Valiant 2 -70 cm

Versatile adjustability



Partnumber: 938828

Entrance plate for Valiant 2 - 70 cm

Even easier entrance to the treadmill



Partnumber: 938829

Handrails, Side - Fixed for extended Valiant 2

Extra support for the test subject



Partnumber: 938820



Robust and reliable treadmill for sports medicine



Specifications

Workload			Connectivity		
Maximum speed	25 km/h	15.5 mph	USB connector	~	
Minimum operational speed	1 km/h	0.6 mph	RS232 in connector	~	
Speed adjustment steps	0.1 km/h	0.1 mph	Dimensions		
Positive elevation	25 %		Screen resolution	800 x 400 pixels	
Elevation adjustment steps	0.5 %		Walking surface length	190 cm	74.8 inch
Optional negative elevation	-10 %		Walking surface width	70 cm	27.6 inch
Accuracy			Step up height	17 cm	6.7 inch
Speed accuracy	5 %		Product length (cm)	255 cm	100.4 inch
Accuracy inclination	0,5 %		Product width (cm)	100 cm	39.4 inch
Comfort			Product height	152 cm	59.8 inch
Allowed user weight	225 kg	496 lbs	Product weight	184 kg	405.7 lbs
User Interface			Power requirements		
English user interface	~		Power cord length	250 cm	98.4 inch
Chinese user interface	~		Power cord IEC 60320 C19 with CEE 7/7 plug	~	
Croatian user interface	~		Power cord NEMA	×	
Czech user interface	~		Power in Horse Power	2.5 HP	
Danish user interface	~		Maximum rated power input	2500 VA	
Dutch user interface	~		115 V AC 50/60 Hz (2 phases)	~	
Finnish user interface	~		230 V AC 50/60 Hz	~	
French user interface	~		Maximum motor power	1.8 kW	
German user interface	~		Standards & Safety		
Greek user interface	~		IEC 60601-1:2005	~	
Hungarian user interface	~		ISO 13485:2016 compliant	~	
Italian user interface	~		ISO 9001:2015 compliant	~	
Japanese user interface	~		Standard emergency lanyard	~	
Korean user interface	~		Certification		
Latvian user interface	~		CE class Im according to MDD93/42/EEC	~	
Lithuanian user interface	~		CE class of product with optional SpO2	lla	
Norwegian user interface	~		CE class of product with optional BPM	lla	
Polish user interface	~		CB according to IECEE CB	~	
Portugese user interface	~				
Romanian user interface	~				
Russian user interface	~				
Spanish user interface	~				
Swedish user interface	~				
Turkish user interface	~				

Order info

Partnumber: 938902

Ukrainian user interface



Distributed by Zernikepark 16 9747 AN Gronin

Pontou 3 11528 Athens Afghanistan Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel: +31 50 5712811
Fax: +31 50 5716746
E-mail: ask@lode.nl
Internet: www.lode.nl

^{*}Specifications are subject to change without notice.